

## Lower Extremity Physiologic Study Single Level (Vantage ABI System Exam)

XYZ Diagnostics, Inc.  
225 Wireless Blvd.  
Anywhere, NY 11788  
1-888-95-TESTS

Patient Name: Mr John Smith    Patient Sex: M    Patient ID: 1234567890    Age: 70    Exam Date: 11/09/10 14:52:29  
Patient DOB: 03/10/1940

Risk Factors	
<input checked="" type="checkbox"/> Tobacco Use	<input type="checkbox"/> Hypertension
<input type="checkbox"/> Diabetes	<input type="checkbox"/> Hyperlipidemia
<input checked="" type="checkbox"/> Heart Disease	<input type="checkbox"/> Stroke/TIA
<input type="checkbox"/> Previous CV Event	<input type="checkbox"/> Prev. Vascular Surgery

Current Symptoms											
<input checked="" type="checkbox"/> Intermittent Claudication	<table border="1"> <thead> <tr> <th colspan="2">Location</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/> Buttocks</td> <td></td> </tr> <tr> <td><input type="checkbox"/> Thigh</td> <td></td> </tr> <tr> <td><input type="checkbox"/> Calf</td> <td></td> </tr> <tr> <td><input checked="" type="checkbox"/> Feet</td> <td></td> </tr> </tbody> </table>	Location		<input type="checkbox"/> Buttocks		<input type="checkbox"/> Thigh		<input type="checkbox"/> Calf		<input checked="" type="checkbox"/> Feet	
Location											
<input type="checkbox"/> Buttocks											
<input type="checkbox"/> Thigh											
<input type="checkbox"/> Calf											
<input checked="" type="checkbox"/> Feet											
<input type="checkbox"/> Numbness, Tingling in Feet											
<input type="checkbox"/> Ulcerations											
<input type="checkbox"/> Rest Pain											
<input type="checkbox"/> Gangrene											

ABI Severity of Disease	
> 1.4	Non-compressible
1.0 - 1.4	Normal
0.91 - 0.99	Borderline
0.41 - 0.90	Moderate to Mild
0.00 - 0.40	Severe

ACCF/AHA Guidelines for Management of patients with P.A.D., 2011

Interpreting Physician: Dr. John M. Smithern

Referring Physician: Dr. John M. Smithern

Comments: No evidence of P.A.D.

**RIGHT**

**LEFT**

Brachial

Brachial

**120**  
mmHg

**125**  
mmHg

Ankle


Ankle

**141**  
mmHg

**139**  
mmHg


**1.12**  
Index

**1.11**  
Index




PVR X1

11/09/10 14:51:55



PVR X1

11/09/10 14:51:34



Interpretation: Normal ABI Results

**WALLACH**<sup>®</sup>  
SURGICAL DEVICES  
**SummitDoppler**<sup>™</sup>  
a WALLACH brand

Wallach Surgical Devices  
95 Corporate Drive • Trumbull, CT 06611 USA  
Phone: 800-243-2463 • (203) 799-2000  
Fax: (203) 799-2002  
customerservice@wallachsurgical.com  
www.wallachsurgical.com