

Marlow Nu-Tip® Laparoscopic Instruments

Combining the convenience of disposables
with the cost savings of reusables.

- The quality and feel of a reusable instrument
- Consistent standard of care for every patient
- Ultimate flexibility
- Cost effective

Marlow Nu-Tip® Laparoscopic Instruments

Combining the convenience of disposables with the cost savings of reusables.

The concept is simple. Reuse the handle, change the tip.

The result. A fresh quality instrument for every procedure.

This cost effective approach gives you ultimate flexibility throughout the procedure.

The Quality and Feel of a Reusable Instrument

These instruments provide quality craftsmanship, strength and precise performance. Made of the highest quality stainless steel, the reusable handle/shaft is designed and engineered to deliver the desired feel of a totally reusable instrument.

Consistent Standard of Care for Every Patient

The Disposable Nu-Tip® screws easily onto the reusable handle/shaft. A full array of graspers and dissectors enable you to create your own instrument, while the extremely sharp double acting scissor jaws provide a sharp cut every time. The result is a precision instrument that provides consistent performance for each and every patient.

Ultimate Flexibility

The Nu-Tip® System, featuring disposable stainless steel tips and a reusable handle/shaft, provides you with the opportunity to create your own scissor, grasper or dissector instrument as required. Re-use the handle, change the tip.

Cost Effective

Reusable scissors seldom remain sharp and fully disposable scissors can be costly and add to the waste stream. The Nu-Tip® System offers a hybrid system resulting in instruments that cut costs and minimize instrument inventory.

ORDERING INFORMATION

Reusable Handle/ Shaft Assembly

Item #M2100M	Insulated, Monopolar, 32cm
Item #M2100RTI	Insulated Shaft, Ratchet, 32cm
Item #M2200M	Insulated, Monopolar, 45cm

Single-Use Disposable Tips

Item #M2101D	Scissors, Curved Metzenbaum, Box of 12
Item #M2102D	Scissors, Straight, Box of 12
Item #M2103D	Scissors, Hook, Box of 12
Item #M3104D	Grasper, Alligator, Box of 12
Item #M3111D	Grasper, Babcock, Box of 12
Item #M3112D	Grasper, Pollack (Atraumatic), Box of 12
Item #M5101D	Dissector, Maryland, Box of 12
Item #M5103D	Dissector, Dolphin Nose, Box of 12

CooperSurgical

95 Corporate Drive, Trumbull, CT 06611
Phone: 800.243.2974 • Fax: 800.262.0105
International Phone: 203.601.9818 • Fax: 203.601.4747 • www.coopersurgical.com
©2006 CooperSurgical, Inc. All rights reserved.

US Patent No. 5,368,606 & No.5,571,137 • Nu-Tip® is a registered trademark of CooperSurgical, Inc.
Form # 80098 01/06



EC REP

Leisegang Feinmechanik GmbH
Leibnizstrasse 32
D-10625, Berlin GERMANY