

INCA[®]

Infant Nasal CPAP Assembly

Your Connection for Neonatal CPAP and Bubble CPAP Systems

Unique Benefits

- Five anatomically proportional nasal cannula sizes for a custom fit
- Clear, extra soft cannula helps maintain tissue integrity
- Stay-flex tubing provides custom shaping and maximum stability
- Foam block retainers provide secure placement and stability
- Chin strap accessory minimizes mouth leaks
- Available in both complete and replacement Sets
- Accessory Pack provides extra knit cap sizes, chin straps, hook fasteners, and foam blocks



“Early bubble CPAP reduced delivery room intubations, (the number of) days on mechanical ventilation, postnatal steroid use and was associated with increased postnatal weight gain without any increased complications.”

- Journal of Perinatology 2003 (23):195-199

INCA® Infant Nasal CPAP Assembly

INCA Complete Set

consists of a nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters, one in-line "T" pressure adapter, a knit cap with hook/loop strap, two foam blocks, and a hook/loop chin strap. One INCA Sizing Gauge is included in every box of five sets.



INCA Replacement Set

consists of a nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters, and one in-line pressure "T" adapter.



INCA Accessory Pack

includes five knit caps with hook/loop straps, five chinstraps with ten hook fasteners (not attached), and ten foam blocks.




Nasal Sizing Gauge

For appropriate nasal prong fitting, use the INCA Sizing Gauge. Wipe gauge with sterile 70% isopropyl alcohol prior to each use. Hold the gauge up to the nasal area to determine appropriate cannula size. **Note: Do not place sizing gauge in neonate's nares.**



Simple Bubble CPAP Setup and Components



INCA® Infant Nasal CPAP Assembly—your connection for Neonatal CPAP (Continuous Positive Air Pressure) and Bubble CPAP Systems.

INCA Infant Nasal CPAP Assembly Ordering Information

Product Code	Description
INCA Complete Set	
44-2707	Complete Set with 7.5 Fr Cannula (5 per Box)
44-2709	Complete Set with 9 Fr Cannula (5 per Box)
44-2710	Complete Set with 10.5 Fr Cannula (5 per Box)
44-2712	Complete Set with 12 Fr Cannula (5 per Box)
44-2715	Complete Set with 15 Fr Cannula (5 per Box)
INCA Replacement Set	
44-0707	Replacement Set with 7.5 Fr Cannula (5 per Box)
44-0709	Replacement Set with 9 Fr Cannula (5 per Box)
44-0710	Replacement Set with 10.5 Fr Cannula (5 per Box)
44-0712	Replacement Set with 12 Fr Cannula (5 per Box)
44-0715	Replacement Set with 15 Fr Cannula (5 per Box)
INCA Accessory Pack	
44-2001	Accessory Pack, Extra-Small (5 XS Caps, 10 Foam Blocks)
44-2002	Accessory Pack, Small (5 S Caps, 10 Foam Blocks)
44-2003	Accessory Pack, Large (5 L Caps, 10 Foam Blocks)
44-2004	Accessory Pack, Extra-Large (5 XL Caps, 10 Foam Blocks)
INCA Nasal Sizing Gauge	
025-001	Nasal Sizing Gauge (1 per Pack)

CooperSurgical

CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products, which aid clinicians in the management and treatment of commonly seen conditions.

800.243.2974 | 203.601.5200 | coopersurgical.com | 95 Corporate Drive | Trumbull, CT 06611
Also on YouTube.com/coopersurgical