



**TAMPA CATHETER**

The simple, cost-effective choice for ultrasound examination of the uterus

**Features**

- Unique introducer sheath facilitates insertion into cervical os with minimal trauma.
- 7cm marker provides optimal placement of catheter tip in uterus.
- Balloonless configuration is ideal for post-menopausal patients with stenotic cervical canals, as well as for many nulliparous patients.
- Small catheter size allows easy placement next to ultrasound probe.



**CooperSurgical**

*Keeping You at the Forefront of Women's Health Care™*


# TAMPA CATHETER

The simple, cost-effective choice for  
ultrasound examination of the uterus

## Ordering information

### Tampa Catheter

- PRODUCT NO. 61-2005  
5 French, clear, flexible catheter with  
introducer sheath, 33cm long.
- PACKAGING  
Sterile, 20 sets per box  
4 boxes (80 sets) per case



## Stainless Steel Instruments for Gynecologic Procedures

### Vu-More™ Wide Rear-Opening Speculum


- Permits easy removal from the patient without disturbing the catheter or instrument
- The longer blade makes it suitable for obese patients

#### Product No. 63-902

Medium

#### Product No. 63-900

Large



### Standard Luer Needle Extension


- Allows use of less expensive standard size needles for paracervical blocks
- Material can be bent to angle needed

#### Product No. 433-131A

Reusable 3" extension

#### Product No. 433-132

Reusable 4" extension



### Weisman-Graves Open-Sided Speculum

- Speculum design permits easy removal without disturbing catheter or instruments
- Ideal time-saving device for Hysterosonography
- Easily adjustable to most anatomic needs
- Available with either left or right side opening

#### Product No. 391-667A

Left side open – Medium

#### Product No. 391-668A

Right side open – Medium

#### Product No. 391-667

Left side open – Large

#### Product No. 391-668

Right side open – Large



### Packaging

All of the products listed above are packaged individually in a clear plastic bag. Four per box. Non-sterile.

TAMPA CATHETER



CooperSurgical

Keeping You at the Forefront  
of Women's Health Care™