

Belly Bands™

Abdominal Fetal Monitoring Belt

Features

- **Flexible and easy to use** — allows mother to place band comfortably over abdomen
- **Comes in pre-cut versions** or in a 10 yard roll
- **14 inch pre-cut length** — covers the belly and secures fetal monitor placement
- **Breathable material** — helps wick away moisture




Belly Bands™

Abdominal Fetal Monitoring Belt

Instructions

- Choose appropriate size
- Open Belly Band and gather with both hands
- Slip Belly Band over both feet
- Slip Belly Band comfortably over abdomen



ORDERING INFORMATION

13010	10" Width Belly Bands on Roll (10" x 10yds.)	1 Roll
13011	10" Width Belly Bands - Pre-cut (10" x 14")	25 per Box
13012	12" Width Belly Bands on Roll (12" x 10yds.)	1 Roll
13013	12" Width Belly Bands - Pre-cut (12" x 14")	25 per Box
13014	14" Width Belly Bands on Roll (14" x 10yds.)	1 Roll
13015	14" Width Belly Bands - Pre-cut (14" x 14")	25 per Box
13016	16" Width Belly Bands on Roll (16" x 10yds.)	1 Roll
13017	16" Width Belly Bands - Pre-cut (16" x 14")	25 per Box
13019	18" Width Belly Bands - Pre-cut (18" x 14")	25 per Box



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products which aid clinicians in the management and treatment of commonly seen conditions.

203.601.5200 | 800.243.2974 | coopersurgical.com | 95 Corporate Dr. | Trumbull, CT 06611
Also on YouTube at www.YouTube.com/CooperSurgical

© 2011 CooperSurgical, Inc.

Form #81238 • 11/11