

ALLY Uterine Positioning System™

Facilitating Static Control of Anatomy
During Gynecologic Laparoscopic Procedures



Provides Constant Static Control and Greater Exposure Where and When it's Needed

The ALLY UPS system provides exceptional static control, creating greater exposure and enhancing visualization during all GYN laparoscopic procedures. From Myomectomy, to TLH and even Sacrocolpopexy, having the ALLY UPS manipulator increases efficiencies while eliminating manpower fatigue.

ALLY UPS ensures steady cephalad pressure, helping to increase the distance between the ureters and other critical anatomy during dissection. Establishing these landmarks and distancing vital anatomy can prevent complications while performing the colpotomy incision.



Hoyte Sacro Tips for Vaginal Prolapse Repair Are Ideally Suited for Use With the ALLY UPS

Adapter and drape are integrated into one device to make draping simple and fast

Adapter can be used with either the RUMI II, Advincula Arch and Advincula Delineator Uterine Manipulators



ALLY UPS console connects efficiently to the right side of the OR table using a standardized bed rail clamp

Foot pedal operation provides quick response and accurate positioning.

ALLY Uterine Positioning System™

Ensures Hands-Free Static Control of Tissue During All GYN Laparoscopic Procedures



The ALLY Uterine Positioning System (UPS) creates exposure and static control of the anatomy during operative and diagnostic GYN laparoscopy. This table-mounted system is designed to provide a full range of motion.

The ALLY UPS is designed to mount, position and secure RUMI® II, and Advincula Arch™ and Advincula Delineator uterine manipulators during laparoscopic procedures. Features and benefits include:

- Securely holds uterus in position for optimal visualization, access and patient safety
- Frees your assistant to perform other tasks during the procedure
- Easy-to-use foot pedal release facilitates positioning
- Connects to any operating table
- Fits easily within a limited workspace, offering easy access with greater efficiency





ALLY Uterine Positioning System Cart

- This unique cart is designed to bring the ALLY UPS directly to the OR table
- Offers safe and reliable maneuverability throughout the hospital
- Provides fast set up and removal of the ALLY UPS

Advincula Delineator™
with Koh-Efficient® Technology



Advincula Arch
with Uterine Tip and Koh-Efficient®



RUMI II
with Uterine Tip and Koh-Efficient



Set Up Easy As 1, 2, 3...

CONNECT



Once the patient is positioned and legs are up in the stirrups, the ALLY UPS is mounted behind the stirrup bracket.

DRAPE



Attach the proper adapter drape to the ALLY UPS and slide the drape down the neck of the ALLY UPS.

DOCK



Insert the uterine manipulator into the patient. Align the adapter to the uterine manipulator and attach.



ALLY UPS



ALLY UPS Cart
(shown with ALLY UPS)



ALLY UPS
Adapter
Including Drape



Advincula
Delineator
Adapter



Advincula Delineator



RUMI II with Uterine Tip and Koh-Efficient



Advincula Arch
with Uterine Tip and Koh-Efficient

ALLY Uterine Positioning System®

Product Code	Description
AU-UPS	ALLY Uterine Positioning System (Domestic)
AU-INT	ALLY Uterine Positioning System (International)
AU-CART	ALLY Uterine Positioning System Cart
AU-AD	ALLY UPS Adapter Drape for RUMI II and Advincula Arch (20/Box)
UMH650	RUMI II Uterine Manipulator Handle
UMH750	Advincula Arch Uterine Manipulator Handle
AU-AD-DLNTR	ALLY UPS Adapter Drape for Delineator (20/Box)

Advincula Delineator™ with Ultem® Koh-Efficient®

Product Code	Description
AD750-KE25	Advincula Delineator 2.5CM Ultem Plastic cup (3/Box)
AD750-KE30	Advincula Delineator 3.0CM Ultem Plastic cup (3/Box)
AD750-KE35	Advincula Delineator 3.5CM Ultem Plastic cup (3/Box)
AD750-KE40	Advincula Delineator 4.0CM Ultem Plastic cup (3/Box)

Advincula Delineator with Stainless Steel Koh-Efficient®

Product Code	Description
AD750-KES25	Advincula Delineator 2.5CM Stainless Steel cup (3/Box)
AD750-KES30	Advincula Delineator 3.0CM Stainless Steel cup (3/Box)
AD750-KES35	Advincula Delineator 3.5CM Stainless Steel cup (3/Box)
AD750-KES40	Advincula Delineator 4.0CM Stainless Steel cup (3/Box)

Advincula Delineator™ with Flexible Soft Cup Koh-Efficient®

Product Code	Description
AD750SC-KE25	Advincula Delineator 2.5CM Flexible Soft cup (3/Box)
AD750SC-KE30	Advincula Delineator 3.0CM Flexible Soft cup (3/Box)
AD750SC-KE35	Advincula Delineator 3.5CM Flexible Soft cup (3/Box)
AD750SC-KE40	Advincula Delineator 4.0CM Flexible Soft cup (3/Box)