

ER-Spec® | by CooperSurgical®

A BRILLIANT SOLUTION



Cordless high-powered illumination for 30+ minutes that doesn't generate damaging heat.



Save time and increase efficiency – no assembly or extra equipment required.



Single-use speculum reduces cross contamination and eliminates reprocessing costs.



Available in 5 sizes to accommodate varied patient anatomies.

Ordering Information

NUMBER	DESCRIPTION
C020130-1	ER-Spec® (Extra-Small)
C020100-1	ER-Spec® (Small)
C020110-1	ER-Spec® (Medium)
C020120-1	ER-Spec® (Large)
C020140-1	ER-Spec® (Extra-Large)
C020300	ER-Spec® (Small - Sterile)
C020310	ER-Spec® (Medium - Sterile)
C020320	ER-Spec® (Large - Sterile)

Integrated light source



5 sizes available

[Learn More](#)

IMPORTANT SAFETY INFORMATION ER-SPEC® is a **Single-Use** lighted vaginal specula indicated for use on women undergoing procedures requiring vaginal access and exposure. The device uses opposing bi-valve plastic blades to separate the vaginal wall, exposing the interior of the vagina and exterior of the cervix. **Warnings: Never** use with laser equipment, gamma radiation, autoclave sterilization. Refer to Instructions for Use, for detailed instructions and risk information www.coopersurgical.com/lighted-speculum-instructions-for-use

800.243.2974 | www.coopersurgical.com



CooperSurgical®