

Her-Spec® | by CooperSurgical®

A BRILLIANT SOLUTION



Cordless high-powered illumination for 30+ minutes that doesn't generate damaging heat.



Save time and increase efficiency – no assembly or extra equipment required.



Single-use speculum reduces cross contamination and eliminates reprocessing costs.




CLICK-FREE design eliminates sound associated with ratchet mechanisms.



Ordering Information

NUMBER	DESCRIPTION
C020500	Her-Spec® (Small)
C020510	Her-Spec® (Medium)
C020520	Her-Spec® (Large)
C020600	Her-Spec® (Small - Sterile)
C020610	Her-Spec® (Medium - Sterile)
C020620	Her-Spec® (Large - Sterile)



[Learn More](#)

IMPORTANT SAFETY INFORMATION HER-SPEC® is a **Single-Use** lighted vaginal specula indicated for use on women undergoing procedures requiring vaginal access and exposure. The device uses opposing bi-valve plastic blades to separate the vaginal wall, exposing the interior of the vagina and exterior of the cervix. **Warnings: Never** use with laser equipment, gamma radiation, autoclave sterilization. Refer to Instructions for Use, for detailed instructions and risk information www.coopersurgical.com/lighted-speculum-instructions-for-use

800.243.2974 | www.coopersurgical.com



CooperSurgical®