

# A BRILLIANT SOLUTION



Cordless high-powered illumination for 30+ minutes that doesn't generate damaging heat.



Save time and increase efficiency – no assembly or extra equipment required.



Single-use speculum reduces cross contamination and eliminates reprocessing costs.




Side-opening design is engineered for easier instrumentation access.

## Disposable, Bright Built-In Light

A single-use, side-opening vaginal speculum with integrated LED light source.

## Ordering Information

NUMBER	DESCRIPTION
C020001	OfficeSpec® (Side-Opening Medium)



[Learn More](#)

**IMPORTANT SAFETY INFORMATION OFFICESPEC®** is a **Single-Use** lighted vaginal specula indicated for use on women undergoing procedures requiring vaginal access and exposure. The device uses opposing bi-valve plastic blades to separate the vaginal wall, exposing the interior of the vagina and exterior of the cervix. **Warnings: Never** use with laser equipment, gamma radiation, autoclave sterilization. Refer to Instructions for Use, for detailed instructions and risk information [www.coopersurgical.com/lighted-speculum-instructions-for-use](http://www.coopersurgical.com/lighted-speculum-instructions-for-use)

800.243.2974 | [www.coopersurgical.com](http://www.coopersurgical.com)



CooperSurgical®