

# Euro-Med™ Stainless Steel Instruments

Precision engineered  
for superb handling





CooperSurgical®

## Euro-Med™, Excellence in Gynecological Specialty Instruments

---

Designed for enhanced control and precision, Euro-Med™ by CooperSurgical offers a comprehensive line of German-engineered stainless steel specula, biopsy punches, and other specialty instruments for women's healthcare.





SCAN TO  
LEARN MORE

# Euro-Med™ Stainless Steel Biopsy Punches

Blue markings indicate bite shape.



## A. Tischler-Morgan™ Biopsy Punch €€

The standard bite incorporates a single tooth on the upper and lower jaws that draws tissue in for sample retention.

| NUMBER  | LENGTH | BITE  |
|---------|--------|---|
| 64-689  | 21 cm  | 7 mm x 3 mm  |
| 64-689L | 25 cm  | 7 mm x 3 mm  |

## B. Baby Tischler Biopsy Punch €€

Small bite incorporates a tooth permitting sampling at the endocervical canal.

| NUMBER  | LENGTH | BITE  |
|---------|--------|---|
| 64-691  | 21 cm  | 4.1 mm x 2.3 mm  |
| 64-691L | 26 cm  | 4.1 mm x 2.3 mm  |

## C. Tischler-Kevorkian™ Biopsy Punch €€

Large rectangular bite includes upper and lower teeth.

| NUMBER | LENGTH | BITE  |
|--------|--------|---|
| 64-679 | 22 cm  | 9 mm x 3 mm  |

## D. Kevorkian-Pacific™ Biopsy Punch €€

A slightly smaller bite and additional anchoring teeth on the bottom jaw.

| NUMBER | LENGTH | BITE  |
|--------|--------|---|
| 64-675 | 22 cm  | 9 mm x 3 mm  |

## E. Mini-Townsend™ Biopsy Punch €€

| NUMBER | LENGTH | BITE  |
|--------|--------|---|
| 64-690 | 21 cm  | 4 mm x 2.3 mm  |

## F. Coppleson Biopsy Punch €€

The unique jaw shape yields a distinct sample orientation for histologic analysis.


| NUMBER | LENGTH | BITE  |
|--------|--------|---|
| 64-660 | 22 cm  | 7.5 mm x 3 mm  |

## G. Baggish Biopsy Punch €€

This punch provides an oval-shaped bite.

The jaws have two bottom teeth and a serrated bed.

| NUMBER | LENGTH | BITE  |
|--------|--------|---|
| 64-695 | 22 cm  | 6.5 mm x 3.5 mm  |






# Euro-Med™ Stainless Steel Biopsy Punches

## H. Burke™ Biopsy Punch

The design features an open handle and an oval-shaped bite for small samples.

| NUMBER | LENGTH | BITE  |
|--------|--------|---|
| 64-677 | 22 cm  | 5.5 mm x 3 mm  |



## I. Euro-Med™ Keyes Vulvar Punch Set and Singles

This set includes six stainless steel interchangeable tips ranging from 3 mm to 8 mm in diameter, and a universal handle. The individual tips can be placed on the threaded handle as required.

| NUMBER | LENGTH   |
|--------|--|
| 68-180 | Set includes: five tips, storage tray and wooden box |



# Euro-Med™ Rotating Biopsy Punch Handles

## A 360-Degree Rotating Shaft for Sample Collection

These rotating punch tips fit either the Euro-Med™ closed rotating handle or the Townsend open rotating handle.



Euro-Med™ Rotating Handle



Townsend™ Rotating Handle

### Tischler-Morgan™ Tip €€

| NUMBER | LENGTH | BITE        |                                     |
|--------|--------|-------------|-------------------------------------|
| 64-484 | 23 cm  | 7 mm x 3 mm | <input checked="" type="checkbox"/> |

### Mini Tip Down €€

| NUMBER | LENGTH | BITE        |                                     |
|--------|--------|-------------|-------------------------------------|
| 64-656 | 23 cm  | 7 mm x 3 mm | <input checked="" type="checkbox"/> |

### Tischler Tip Up €€

| NUMBER | LENGTH | BITE        |                                     |
|--------|--------|-------------|-------------------------------------|
| 64-751 | 23 cm  | 7 mm x 3 mm | <input checked="" type="checkbox"/> |

### Mini-Townsend™ €€

| NUMBER | LENGTH | BITE          |                                     |
|--------|--------|---------------|-------------------------------------|
| 64-460 | 23 cm  | 4 mm x 2.3 mm | <input checked="" type="checkbox"/> |

### Tischler Tip Down €€

| NUMBER | LENGTH | BITE        |                                     |
|--------|--------|-------------|-------------------------------------|
| 64-756 | 23 cm  | 7 mm x 3 mm | <input checked="" type="checkbox"/> |

### Oval Tip €€

| NUMBER | LENGTH | BITE        |                                     |
|--------|--------|-------------|-------------------------------------|
| 64-489 | 23 cm  | 5 mm x 3 mm | <input checked="" type="checkbox"/> |

### Mini Tip Up €€

| NUMBER | LENGTH | BITE        |                                     |
|--------|--------|-------------|-------------------------------------|
| 64-651 | 23 cm  | 7 mm x 3 mm | <input checked="" type="checkbox"/> |

### Kevorkian Tip €€

| NUMBER | LENGTH | BITE        |                                     |
|--------|--------|-------------|-------------------------------------|
| 64-488 | 23 cm  | 8 mm x 3 mm | <input checked="" type="checkbox"/> |


# Euro-Med™ Standard Biopsy Punches



## Eppendorfer Tip €€

This oval cutting jaw obtains an 8 mm x 4.5 mm cervical specimen. A basket in the jaw helps retain the biopsy sample.


| NUMBER | LENGTH | BITE          |  |
|--------|--------|---------------|--|
| 64-687 | 21 cm  | 8 mm x 4.5 mm |  |



## Wittner Tip €€

The angle of the jaw enables extra access to areas of the cervix. The large teeth and strong bite provide a maximum tissue sample area.

| NUMBER | LENGTH  | BITE        |  |
|--------|---------|-------------|--|
| 64-681 | 21.5 cm | 8 mm x 4 mm |  |



## Kevorkian Tip €€

This punch facilitates simple insertion and sampling of the cervical canal.

| NUMBER | LENGTH  | BITE        |  |
|--------|---------|-------------|--|
| 64-685 | 26.5 cm | 9 mm x 3 mm |  |



## Schubert Tip €€

Features a large bite and an angled handle.


| NUMBER | LENGTH  | BITE          |  |
|--------|---------|---------------|--|
| 64-670 | 21.5 cm | 8 mm x 4.5 mm |  |

# Euro-Med™ Stainless Steel Specula

## A. Graves €€

Available with standard yoke and an open-sided configuration. The Graves blade is also available in longer lengths with a Vu-More™ or Vu-Max™ yoke.

| NUMBER | DESCRIPTION | SIZE MM  |
|--------|-------------|----------|
| 64-100 | Small       | W19 L82  |
| 64-102 | Medium      | W33 L114 |
| 64-104 | Large       | W35 L127 |
| 64-108 | Open-sided  | W35 L100 |




A.

## B. Pederson

Pederson vaginal speculum is a double-bladed instrument with a screw mechanism to keep the blades open once inserted into the vaginal canal. The speculum is designed to aid in the visualization of the vaginal canal for diagnostic and/or procedural purposes.

| NUMBER  | DESCRIPTION  | SIZE MM  |
|---------|--------------|----------|
| 391-671 | Virgin       | W13 L78  |
| 64-110  | Small        | W13 L82  |
| 391-674 | Extra-narrow | W17 L108 |
| 64-112  | Medium       | W22 L114 |
| 64-114  | Large        | W25 L130 |




B.

## C. Vu-More™ €€

The Vu-More provides a 25 percent larger viewing area at the yoke.

| NUMBER | DESCRIPTION | SIZE MM  |
|--------|-------------|----------|
| 64-130 | Mini        | W34 L93  |
| 63-902 | Medium      | W33 L113 |
| 63-900 | Large       | W35 L127 |
| 63-905 | X Large     | W33 L150 |
| 63-906 | 2X Large    | W33 L160 |
| 63-907 | 3X Large    | W33 L170 |
| 63-908 | 4X Large    | W33 L180 |



C.



# Euro-Med™ Stainless Steel Specula



## Vu-Max™ C€

Offered in both Pederson and Graves configurations.

| NUMBER | DESCRIPTION     | SIZE MM  |
|--------|-----------------|----------|
| 64-113 | Pederson Medium | W22 L114 |
| 64-115 | Pederson Large  | W25 L120 |
| 64-218 | Graves Medium   | W33 L114 |
| 64-219 | Graves Large    | W35 L125 |
| 64-225 | X Large         | W35 L150 |
| 64-226 | 2X Large        | W35 L160 |
| 64-227 | 3X Large        | W35 L170 |
| 64-228 | 4X Large        | W35 L180 |



## Snowman Specula™

Longer blades expose access to the cervix.

| NUMBER | DESCRIPTION                | SIZE MM     |
|--------|----------------------------|-------------|
| 64-107 | Graves Blade               | W35 L152 C€ |
| 64-117 | Pederson Blade             | W25 L152 C€ |
| 64-127 | One each/Graves & Pederson |             |



## ExpandableView™ Integrated Speculum and Lateral Retractor

Combines a traditional speculum with a lateral wall retractor.

| NUMBER | DESCRIPTION                       | SIZE MM     |
|--------|-----------------------------------|-------------|
| 64-161 | Medium                            | W22 L109 C€ |
| 64-162 | Large                             | W22 L133 C€ |
| 64-163 | Stainless Steel Instrument Holder |             |



## Auvard C€

This 2.5-pound weighted speculum for posterior vaginal wall retraction in vaginal procedures.

| NUMBER | DESCRIPTION            | SIZE MM |
|--------|------------------------|---------|
| 64-422 | Med (Wide/Short Blade) | W43 L55 |
| 64-421 | Lg (Narrow/Long Blade) | W35 L90 |



## Weisman Graves Side Opening C€

Available in two sizes with left or right openings.

| NUMBER   | DESCRIPTION          | SIZE MM  |
|----------|----------------------|----------|
| 391-667A | Left side open, Med  | W33 L102 |
| 391-667  | Left side open, Lg   | W36 L121 |
| 391-668A | Right side open, Med | W35 L102 |
| 391-668  | Right side open, Lg  | W36 L120 |




## Pederson Side Opening C€

The Pederson open-sided specula features either a left or right side opening.

| NUMBER  | DESCRIPTION          | SIZE MM  |
|---------|----------------------|----------|
| 64-112L | Left side open, Med  | W22 L114 |
| 64-114L | Left side open, Lg   | W25 L120 |
| 64-112R | Right side open, Med | W22 L114 |
| 64-114R | Right side open, Lg  | W25 L120 |

# Euro-Med™ Stainless Steel Specula

Vu-Max™



YOKE  
Standard  
W: 30 mm



YOKE  
Vu-More™  
W: 40 mm



YOKE  
Vu-Max  
W: 70 mm



## Specula Blades and Yoke Dimensions

Use our speculum sizing guide to select the proper tools for your practice. Specula are available in a variety of sizes.

\*Approximate size and shape.  
See pages 7–8 for a full product selection and model numbers.

### Vu-Max™ Yoke – Graves Blades

Medium – W: 33 mm L: 114 mm  
Large – W: 35 mm L: 125 mm

### Vu-Max Yoke – Pederson Blades

Medium – W: 22 mm L: 114 mm  
Large – W: 25 mm L: 120 mm

### C. Vu-More™ Yoke – Graves Blades

Mini – W: 34 mm L: 93 mm  
Medium – W: 33 mm L: 113 mm  
Large – W: 35 mm L: 127 mm

### D. Standard Yoke – Pederson Blades

Virgin – W: 13 mm L: 78 mm  
Small – W: 13 mm L: 82 mm  
Extra Narrow – W: 17 mm L: 108 mm  
Medium – W: 22 mm L: 114 mm  
Large – W: 25 mm L: 130 mm

Also Available with  
Side Opening Yoke

### E. Standard Yoke – Graves Blades

Small – W: 19 mm L: 82 mm  
Medium – W: 33 mm L: 114 mm  
Large – W: 35 mm L: 127 mm

Also Available with  
Side Opening Yoke

### F. Vu-More Yoke – Graves or Pederson Blade

Snowman™ – W: 35 mm L: 152 mm  
Snowman – W: 25 mm L: 152 mm

### G. Vu-Max and Vu-More Yokes – Graves Blades

X Large – W: 35 mm L: 150 mm  
2X Large – W: 35 mm L: 160 mm  
3X Large – W: 35 mm L: 170 mm  
4X Large – W: 35 mm L: 180 mm




# Euro-Med™ Stainless Steel Lateral Wall Retractors



A.

## A. Cer-View™ Lateral Wall Retractor CA★ CE

| NUMBER | DESCRIPTION | SIZE MM  |
|--------|-------------|----------|
| 64-310 | Medium      | W22 L150 |
| 64-320 | Large       | W22 L160 |



B.

## B. Cer-View™ Lateral Wall Retractor CE

| NUMBER | DESCRIPTION | SIZE MM  |
|--------|-------------|----------|
| 64-307 | Extra Large | W32 L190 |



C.

## C. Jackson Vaginal Retractor CE

| NUMBER | DESCRIPTION | SIZE MM  |
|--------|-------------|----------|
| 64-322 | Medium      | W40 L98  |
| 64-324 | Large       | W40 L115 |

# Euro-Med™ Stainless Steel Uterine Sounds and Dilators

## A. Dil-Os™ Uterine Sound/Dilator

| NUMBER | DESCRIPTION                 | SIZE          |
|--------|-----------------------------|---------------|
| 64-336 | Slight Curve                | L12.7"/322 mm |
| 64-337 | Strong Curve                | L12.5"/317 mm |
| 64-338 | Set of two (64-336, 64-337) |               |



## B. Sims Uterine Sound

| NUMBER | DESCRIPTION                             | SIZE       |
|--------|---|------------|
| 64-600 | Sims Uterine Sound Silver-Plated Copper | L12.5"/317 |
| 64-601 | Sims Uterine Sound Stainless Steel      | L12.5"/317 |

# Euro-Med™ Stainless Steel Uterine Sounds and Dilators



## Pratt Dilator Set €€

Available individually or in a set of 8.

| NUMBER  | SIZE                              |
|---------|-----------------------------------|
| 391-830 | 13 Fr/15 Fr – 41 Fr/43 Fr (8/Set) |
| 391-833 | 13 Fr/15 Fr                       |
| 391-834 | 17 Fr/19 Fr                       |
| 391-835 | 21 Fr/23 Fr                       |
| 391-836 | 25 Fr/27 Fr                       |
| 391-837 | 29 Fr/31 Fr                       |
| 391-838 | 33 Fr/35 Fr                       |
| 391-839 | 37 Fr/39 Fr                       |
| 391-840 | 41 Fr/43 Fr                       |
| 391-841 | 45 Fr/47 Fr                       |
| 391-842 | 49 Fr/51 Fr                       |
| 391-843 | 53 Fr/55 Fr                       |
| 391-844 | 57 Fr/59 Fr                       |



## Hegar Dilator Set €€

Stainless steel Hegar dilators provide the option of double-ended dilators. Available individually or in a set of 8.

| NUMBER | SIZE MM  |
|--------|--|
| 61207  | 9 Fr (3 mm)/12 Fr (4 mm)/<br>51 Fr (17 mm)/54 Fr (18 mm) (8/Set) |
| 61208  | 9 Fr (3 mm)/12 Fr (4 mm)   |
| 61209  | 15 Fr (5 mm)/18 Fr (6 mm)  |
| 61210  | 21 Fr (7 mm)/24 Fr (8 mm)  |
| 61211  | 27 Fr (9 mm)/30 Fr (10 mm)                                       |
| 61212  | 33 Fr (11 mm)/36 Fr (12 mm)                                      |
| 61213  | 39 Fr (13 mm)/42 Fr (14 mm)                                      |
| 61214  | 45 Fr (15 mm)/48 Fr (16 mm)                                      |
| 61215  | 51 Fr (17 mm)/54 Fr (18 mm)                                      |



## Hank Cervical Dilators €€

Hank cervical dilators available individually or in a set of six ranging from 9 French (3 mm) to 20 French (6.67 mm). The set comes in a custom-fitted and numbered pouch.

| NUMBER | SIZE MM                         |
|--------|---------------------------------|
| 64-370 | 6 Dilators and Storage Pouch    |
| 61201  | 9 Fr (3 mm)/10 Fr (3.33 mm)     |
| 61202  | 11 Fr (3.67 mm)/12 Fr (4 mm)    |
| 61203  | 13 Fr (4.33 mm)/14 Fr (4.67 mm) |
| 61204  | 15 Fr (5 mm)/16 Fr (5.33 mm)    |
| 61205  | 17 Fr (5.67 mm)/18 Fr (6 mm)    |
| 61206  | 19 Fr (6.33 mm)/20 Fr (6.67 mm) |



## Mini Cervical Dilator Set €€

| NUMBER | SIZE MM                              |
|--------|--------------------------------------|
| 64-374 | 5 Dilators, Handle and Storage Pouch |



## Small Cervical Dilator Set €€

| NUMBER | SIZE MM                      |
|--------|------------------------------|
| 64-372 | 6 Dilators and Storage Pouch |



# Euro-Med™ Stainless Steel Endocervical Curettes



Kevorkian with Basket

| NUMBER | LENGTH     |
|--------|------------|
| 66-442 | 12"/305 mm |



Kevorkian without Basket

| NUMBER | LENGTH     |
|--------|------------|
| 66-440 | 12"/305 mm |



Townsend™ C€

| NUMBER | LENGTH       |
|--------|--------------|
| 66-445 | 12.3"/312 mm |



Kevorkian Collector™

| NUMBER | LENGTH     |
|--------|------------|
| 66-441 | 12"/305 mm |



London™

| NUMBER | LENGTH       |
|--------|--------------|
| 64-490 | 11.7"/297 mm |



Endocrette™ C€

| NUMBER | LENGTH       |
|--------|--------------|
| 66-447 | 12.3"/312 mm |

# Euro-Med™ Endometrial Curettes



Handle Style

A.



B.



## A. Novak €€

| NUMBER | SIZE | LENGTH      |
|--------|------|-------------|
| 64-624 | 1 mm | 9.8"/248 mm |
| 64-626 | 2 mm | 9.8"/248 mm |
| 64-628 | 3 mm | 9.8"/248 mm |
| 64-630 | 4 mm | 9.8"/248 mm |

## B. Randall €€

| NUMBER | SIZE   | LENGTH      |
|--------|--------|-------------|
| 64-620 | 3.9 mm | 9.9"/252 mm |



Handle Style

C.

## C. Sims

| NUMBER  | SIZE    | LENGTH        |
|---------|---------|---------------|
| 392-110 | 6 mm    | 11"/280 mm €€ |
| 392-111 | 7 mm    | 11"/280 mm €€ |
| 392-112 | 8 mm    | 11"/280 mm €€ |
| 392-113 | 9 mm    | 11"/280 mm €€ |
| 392-114 | 11.5 mm | 11"/280 mm €€ |
| 392-116 | 17 mm   | 11"/280 mm    |

# Euro-Med™ Stainless Steel Hooks



## Angle Hooks

| NUMBER | DESCRIPTION   | LENGTH      |
|--------|---------------|-------------|
| 36-630 | Straight Hook | 9.8"/249 mm |
| 36-632 | Right Hook    | 9.6"/244 mm |



## Iris Hooks

| NUMBER | LENGTH       |
|--------|--------------|
| 36-625 | 10.3"/262 mm |



## Marlow Corkscrew Hook €€

| NUMBER | LENGTH      |
|--------|-------------|
| 67-451 | 9.3"/236 mm |



## Double IUD Hook

| NUMBER  | DESCRIPTION  | LENGTH      |
|---------|--------------|-------------|
| 392-151 | 10.6"/269 mm | 9.8"/249 mm |



## Emmett Tenaculum Hook €€

| NUMBER | LENGTH      |
|--------|-------------|
| 66-450 | 9.8"/249 mm |



## IUD Hook €€

| NUMBER  | LENGTH       |
|---------|--------------|
| 392-148 | 10.3"/262 mm |



# Euro-Med™ Stainless Steel Cervical Tenacula



**Braun Tenaculum  
Forceps** €€

| NUMBER | LENGTH      |
|--------|-------------|
| 64-420 | 9.9"/252 mm |



**Schroeder  
Tenaculum** €€

| NUMBER | LENGTH      |
|--------|-------------|
| 64-410 | 9.9"/252 mm |



**Duplay Tenaculum  
Double Curved** €€

| NUMBER | LENGTH     |
|--------|------------|
| 61352  | 11"/280 mm |



**Goldstein Grasp  
Cervical Stabilizer** €€

| NUMBER  | LENGTH      |
|---------|-------------|
| 64-6026 | 9.4"/239 mm |

# Euro-Med™ Stainless Steel Dressing and Tissue Forceps



A.

## A. Allis

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61981  | 4 x 5 Teeth | 6"     |
| 61983  | 5 x 6 Teeth | 7.5"   |
| 61882  | 5 x 6 Teeth | 10"    |



B.

## B. Dressing Serrated

| NUMBER | LENGTH |
|--------|--------|
| 61804  | 5.6"   |
| 61805  | 6.2"   |



C.

## C. Potts-Smith Tissue

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61871  | 1 x 2 Teeth | 7"     |
| 61872  | 1 x 2 Teeth | 8.2"   |

# Euro-Med™ Stainless Steel Dressing and Tissue Forceps



**Potts-Smith Dressing** €€

| NUMBER | LENGTH |
|--------|--------|
| 10-202 | 8.25"  |



**Tissue 1 x 2** €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61841  | 1 x 2 Teeth | 5.75"  |
| 61843  | 1 x 2 Teeth | 7"     |
| 61844  | 1 x 2 Teeth | 8"     |



**Adson Dressing Serrated** €€

| NUMBER | DESCRIPTION   | LENGTH |
|--------|---------------|--------|
| 61830  | Serrated Tips | 4.75"  |



**Adson Tissue 1 x 2** €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61863  | 1 x 2 Teeth | 4.75"  |



**Adson Tissue 1 x 2 Serrated** €€

| NUMBER | DESCRIPTION         | LENGTH |
|--------|---------------------|--------|
| 61864  | 1 x 2 Serrated Tips | 4.75"  |



**Adson-Brown Grasping** €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61866  | 7 x 7 Teeth | 4.75"  |

# Euro-Med™ Stainless Steel Endocervical Specula



## Ratchet Handle Speculum ☾☾

Stainless steel endocervical speculum with quick-release lock.

| NUMBER | TIP SIZE    | LENGTH       |
|--------|-------------|--------------|
| 66-437 | Small 3 mm  | 10.6"/270 mm |
| 66-435 | Medium 4 mm | 10.6"/270 mm |
| 66-436 | Large 5 mm  | 10.6"/270 mm |



## Screw/Knob Handle Speculum ☾☾

The screw/knob handle stainless steel speculum adjusts and maintains position using the threaded thumb screw.

| NUMBER | TIP SIZE    | LENGTH       |
|--------|-------------|--------------|
| 66-431 | Small 3 mm  | 10.6"/270 mm |
| 66-430 | Medium 4 mm | 10.6"/270 mm |
| 66-432 | Large 5 mm  | 10.6"/270 mm |



## No-Lock Handle ☾☾

This stainless steel endocervical speculum offers single-handed action. The No-Lock Handle allows for adjustment as necessary.

| NUMBER | TIP SIZE    | LENGTH       |
|--------|-------------|--------------|
| 66-433 | Small 3 mm  | 10.6"/270 mm |
| 67-435 | Medium 4 mm | 10.6"/270 mm |
| 66-434 | Large 5 mm  | 10.6"/270 mm |



# Each Endocervical Speculum Available in Small, Medium or Large Tip

3 mm  
Small Size "Thin Line" Tip



4 mm  
Medium Size Euro-Med™ Tip



5 mm  
Large Size Euro-Med™ Tip



## Foerster Sponge Forceps €€

Featuring oval tips with serrations, the Foerster is used to hold sponges during the procedure.

| NUMBER | LENGTH |
|--------|--------|
| 16-182 | 9.7"   |

# Euro-Med™ Stainless Steel Uterine Forceps



**Bozeman  
Uterine Forceps** €€

| NUMBER  | LENGTH       |
|---------|--------------|
| 392-660 | 10.2"/257 mm |



**Kelly  
Uterine Forceps** €€

| NUMBER  | LENGTH       |
|---------|--------------|
| 417-131 | 12.5"/318 mm |



**Foerster  
Uterine Polyp Forceps** €€

| NUMBER  | DESCRIPTION | LENGTH     |
|---------|-------------|------------|
| 415-140 | Straight    | 10"/254 mm |
| 415-141 | Curved      | 10"/254 mm |

# Euro-Med™ Stainless Steel Needle Holders



**Mayo-Hegar (Wider Jaw)** €€

| NUMBER | LENGTH |
|--------|--------|
| 62809  | 5.5"   |
| 62810  | 6.3"   |
| 62815  | 7.1"   |
| 62820  | 8.1"   |



**Heaney (Curved)** €€

| NUMBER | LENGTH |
|--------|--------|
| 62825  | 8.3"   |



**Crile-Wood (Narrow Jaw)** €€

| NUMBER | LENGTH |
|--------|--------|
| 62840  | 7.1"   |

# Euro-Med™ Stainless Steel Scissors



## Operating Scissors (Straight) €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 62002  | B/B         | 5.5"   |
| 62008  | S/B         | 5.5"   |
| 62014  | S/S         | 5.5"   |
| 62004  | B/B         | 6.5"   |
| 62005  | B/B         | 7.5"   |
| 62010  | S/B         | 6.5"   |
| 62011  | S/B         | 7.23"  |



## Spencer Suture Scissors €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 62164  | Straight    | 3.5"   |



## Metzenbaum Scissors (Straight or Curved) €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 62116  | Straight    | 7"     |
| 62127  | Straight    | 8"     |
| 62120  | Straight    | 9.1"   |
| 62122  | Straight    | 11.2"  |
| 62131  | Curved      | 5.9"   |
| 62118  | Curved      | 7"     |
| 62119  | Curved      | 8.1"   |
| 62124  | Curved      | 9"     |
| 62126  | Curved      | 11"    |



## Mayo Scissors (Straight or Curved) €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 62063  | Straight    | 5.7"   |
| 62065  | Straight    | 6.8"   |
| 62066  | Straight    | 9"     |
| 62067  | Curved      | 5.5"   |
| 62069  | Curved      | 6.75"  |
| 62070  | Curved      | 9"     |



## Iris Scissors (Straight or Curved) €€

| NUMBER | DESCRIPTION  | LENGTH |
|--------|--------------|--------|
| 62140  | Straight S/S | 4.5"   |
| 62143  | Curved S/S   | 4.5"   |

## Committed to Our Customers

At CooperSurgical, we strive for the highest level of expertise and knowledge of our products and to develop an ongoing relationship with each of our customers that is built on trust, quality and consistency. We look to provide service that exceeds not only our customers' expectations, but industry standards as well. CooperSurgical's Customer Service Department provides full-service order entry with extended hours to facilitate order placement. As an added convenience we also offer online order entry. Call now to have your Customer Service Representative set up your account today. Once you are registered for online ordering, you will be able to place orders, track shipments, review prior purchase history and schedule orders for future delivery.

CE Available for sale where CE mark is required

CA Available for sale in Canada

CE and CA labeling is subject to change. Contact Customer Service or your local sales representative for the most up-to-date information.

## 5 Ways to Order

### Local Sales Representative

[www.coopersurgical.com/pages/salesrepsearch.aspx](http://www.coopersurgical.com/pages/salesrepsearch.aspx)

### Email

U.S.: [Orders@coopersurgical.com](mailto:Orders@coopersurgical.com)

Int'l: [InternationalCustomerService@coopersurgical.com](mailto:InternationalCustomerService@coopersurgical.com)

### E-Ordering (U.S. Only)

Call 800.243.2974 to set up your account today.

### Phone Orders

800.243.2974 or 203.601.5200

(8:00 AM – 8:30 PM EST)

Int'l: +1 203.601.9818

### Online

[ecommerce.coopersurgical.com](http://ecommerce.coopersurgical.com)



SCAN TO  
LEARN MORE

## Warranty Information

CooperSurgical, Inc. ("CooperSurgical") warrants that its products offered in this catalog, **other than** single-use products and consumables, will be free from defects in materials and workmanship for one year from the date of purchase.

CooperSurgical warrants that its single-use products and consumables offered in this catalog will, through the earlier of 90 days from the date of purchase of each such product and the expiration date of that product stated on that product's package, be free from defects in materials and workmanship that are discovered prior to or during the single use of that product.

If CooperSurgical determines that a product offered in this catalog fails to conform to the CooperSurgical warranty applicable to that product, CooperSurgical will, as the sole remedy for that failure to conform, repair or replace that product, at CooperSurgical's discretion, free of charge.

To return a product to CooperSurgical, a customer must comply with CooperSurgical's Returned Goods Policy described in this catalog. A customer will not have any remedy if a CooperSurgical product does not conform to the warranty for that product unless that product is returned to CooperSurgical in accordance with that Returned Goods Policy. CooperSurgical will ship returned products that it repaired or replaced under warranty to the customer who returned those products, at CooperSurgical's expense F.O.B. the customer's facility. Under all other circumstances, CooperSurgical will ship returned products to the customer who returned those products at the customer's expense F.O.B. CooperSurgical's facility.

CooperSurgical's warranties, including those in this catalog, do not cover damage caused by misuse, improper care, improper use of chemicals or cleaning methods, loss, theft, servicing or sharpening by non-authorized personnel or negligent or intentional conduct on the part of the owner or user of the product. CooperSurgical's warranties do not cover normal wear and tear or general maintenance, including cleaning and sharpening services required as a result of normal use. Any modifications or changes to a product will void that product's warranty.

THE ABOVE WARRANTIES ARE IN LIEU OF, AND COOPERSURGICAL HEREBY DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, WRITTEN OR ORAL, WITH RESPECT TO COOPERSURGICAL'S PRODUCTS, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO TERMS, CONDITIONS, UNDERSTANDINGS OR AGREEMENTS THAT PURPORT TO MODIFY THE ABOVE WARRANTIES OR THAT MAKE ANY ADDITIONAL WARRANTIES FOR ANY COOPERSURGICAL PRODUCT SHALL HAVE ANY LEGAL EFFECT UNLESS MADE IN WRITING AND SIGNED BY AN AUTHORIZED COOPERSURGICAL CORPORATE OFFICER.

COOPERSURGICAL SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR LOST PROFITS, DAMAGES FROM LOSS OF USE OR LOST DATA, OR INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ITS WARRANTIES OR OTHERWISE FOR ANY CLAIM RELATED TO COOPERSURGICAL'S PRODUCTS, EVEN IF COOPERSURGICAL HAS BEEN ADVISED, KNEW OR SHOULD HAVE KNOWN OF THE POSSIBILITY OF SUCH DAMAGES. COOPERSURGICAL'S LIABILITY WITH RESPECT TO A PRODUCT COVERED BY A WARRANTY OR OTHERWISE SHALL BE LIMITED IN ALL CIRCUMSTANCES TO THE PURCHASE PRICE OF THAT PRODUCT.

## Returned Goods Policy

### 1. Goods will be accepted for return for the following reasons:

- If shipment was made without Customer's authorization or order
- If incorrect goods were shipped
- If defective goods are covered by the standard warranty

### 2. To return goods, you must contact a Customer Service Representative by telephone (1.800.243.2974 or 1.203.601.5200) for a Returned Merchandise Authorization (RMA) number. Items will not be accepted without an RMA number.

#### Please have the following information:

- Reason you wish to return the goods
- Quantity, description, part number, serial number of the goods
- Date of Receipt of Order
- Customer's Purchase Order and the CooperSurgical Invoice Number

**All used instruments must be cleaned and sterilized prior to being shipped for return.**

### 3. Shipment must be sent prepaid by Customer. Freight collect shipments will not be accepted, and goods will be returned to sender.

### 4. If Customer intends to return equipment ordered in error, the following restocking charges and terms will apply:

- There is no restocking fee if goods are returned within 30 days of shipment
- There is a restocking fee of 25% if goods are returned before 31-60 days of shipment
- Goods must be returned unused, in the original carton, and in marketable condition
- Refurbishing and replacement charges will be added to the restocking charges for damaged or missing items
- No return after 60 days
- No refund on sterile, single-use disposable products

**Send All Returns to: CooperSurgical  
95 Corporate Drive  
Trumbull, CT 06611**





CooperSurgical®

800.243.2974  
coopersurgical.com  
95 Corporate Drive, Trumbull, CT 06611