

Euro-Med™ Stainless Steel Instruments

Precision engineered
for superb handling





CooperSurgical®

Euro-Med™, Excellence in Gynecological Specialty Instruments

Designed for enhanced control and precision, Euro-Med™ by CooperSurgical offers a comprehensive line of German-engineered stainless steel specula, biopsy punches, and other specialty instruments for women's healthcare.





SCAN TO
LEARN MORE

Euro-Med™ Stainless Steel Biopsy Punches

Blue markings indicate bite shape.



A. Tischler-Morgan™ Biopsy Punch €€

The standard bite incorporates a single tooth on the upper and lower jaws that draws tissue in for sample retention.

NUMBER	LENGTH	BITE
64-689	21 cm	7 mm x 3 mm 
64-689L	25 cm	7 mm x 3 mm 

B. Baby Tischler Biopsy Punch €€

Small bite incorporates a tooth permitting sampling at the endocervical canal.

NUMBER	LENGTH	BITE
64-691	21 cm	4.1 mm x 2.3 mm 
64-691L	26 cm	4.1 mm x 2.3 mm 

C. Tischler-Kevorkian™ Biopsy Punch €€

Large rectangular bite includes upper and lower teeth.

NUMBER	LENGTH	BITE
64-679	22 cm	9 mm x 3 mm 

D. Kevorkian-Pacific™ Biopsy Punch €€

A slightly smaller bite and additional anchoring teeth on the bottom jaw.

NUMBER	LENGTH	BITE
64-675	22 cm	9 mm x 3 mm 

E. Mini-Townsend™ Biopsy Punch €€

NUMBER	LENGTH	BITE
64-690	21 cm	4 mm x 2.3 mm 

F. Coppleson Biopsy Punch €€

The unique jaw shape yields a distinct sample orientation for histologic analysis.


NUMBER	LENGTH	BITE
64-660	22 cm	7.5 mm x 3 mm 

G. Baggish Biopsy Punch €€

This punch provides an oval-shaped bite.

The jaws have two bottom teeth and a serrated bed.

NUMBER	LENGTH	BITE
64-695	22 cm	6.5 mm x 3.5 mm 




Euro-Med™ Stainless Steel Biopsy Punches

H. Burke™ Biopsy Punch

The design features an open handle and an oval-shaped bite for small samples.

NUMBER	LENGTH	BITE
64-677	22 cm	5.5 mm x 3 mm



I. Euro-Med™ Keyes Vulvar Punch Set and Singles

This set includes six stainless steel interchangeable tips ranging from 3 mm to 8 mm in diameter, and a universal handle. The individual tips can be placed on the threaded handle as required.

NUMBER	LENGTH
68-180	Set includes: five tips, storage tray and wooden box

Euro-Med™ Rotating Biopsy Punch Handles

A 360-Degree Rotating Shaft for Sample Collection

These rotating punch tips fit either the Euro-Med™ closed rotating handle or the Townsend open rotating handle.



Euro-Med™
Rotating Handle



Townsend™
Rotating Handle

Tischler-Morgan™ Tip €€

NUMBER	LENGTH	BITE	
64-484	23 cm	7 mm x 3 mm	<input checked="" type="checkbox"/>

Mini Tip Down €€

NUMBER	LENGTH	BITE	
64-656	23 cm	7 mm x 3 mm	<input checked="" type="checkbox"/>

Tischler Tip Up €€

NUMBER	LENGTH	BITE	
64-751	23 cm	7 mm x 3 mm	<input checked="" type="checkbox"/>

Mini-Townsend™ €€

NUMBER	LENGTH	BITE	
64-460	23 cm	4 mm x 2.3 mm	<input checked="" type="checkbox"/>

Tischler Tip Down €€

NUMBER	LENGTH	BITE	
64-756	23 cm	7 mm x 3 mm	<input checked="" type="checkbox"/>

Oval Tip €€

NUMBER	LENGTH	BITE	
64-489	23 cm	5 mm x 3 mm	<input checked="" type="checkbox"/>

Mini Tip Up €€

NUMBER	LENGTH	BITE	
64-651	23 cm	7 mm x 3 mm	<input checked="" type="checkbox"/>


Kevorkian Tip €€

NUMBER	LENGTH	BITE	
64-488	23 cm	8 mm x 3 mm	<input checked="" type="checkbox"/>

Euro-Med™ Standard Biopsy Punches



Handle for
Eppendorfer, Wittner
and Kevorkian Tips



Handle for
Schubert Tip



Eppendorfer Tip

This oval cutting jaw obtains an 8 mm x 4.5 mm cervical specimen. A basket in the jaw helps retain the biopsy sample.


NUMBER	LENGTH	BITE
64-687	21 cm	8 mm x 4.5 mm



Wittner Tip

The angle of the jaw enables extra access to areas of the cervix. The large teeth and strong bite provide a maximum tissue sample area.

NUMBER	LENGTH	BITE
64-681	21.5 cm	8 mm x 4 mm



Kevorkian Tip

This punch facilitates simple insertion and sampling of the cervical canal.

NUMBER	LENGTH	BITE
64-685	26.5 cm	9 mm x 3 mm



Schubert Tip

Features a large bite and an angled handle.


NUMBER	LENGTH	BITE
64-670	21.5 cm	8 mm x 4.5 mm

Euro-Med™ Stainless Steel Specula

A. Graves €€

Available with standard yoke and an open-sided configuration. The Graves blade is also available in longer lengths with a Vu-More™ or Vu-Max™ yoke.


NUMBER	DESCRIPTION	SIZE MM
64-100	Small	W19 L82
64-102	Medium	W33 L114
64-104	Large	W35 L127
64-108	Open-sided	W35 L100



B. Pederson

Pederson vaginal speculum is a double-bladed instrument with a screw mechanism to keep the blades open once inserted into the vaginal canal. The speculum is designed to aid in the visualization of the vaginal canal for diagnostic and/or procedural purposes.


NUMBER	DESCRIPTION	SIZE MM
391-671	Virgin	W13 L78
64-110	Small	W13 L82
391-674	Extra-narrow	W17 L108
64-112	Medium	W22 L114
64-114	Large	W25 L130



C. Vu-More™ €€

The Vu-More provides a 25 percent larger viewing area at the yoke.

NUMBER	DESCRIPTION	SIZE MM
64-130	Mini	W34 L93
63-902	Medium	W33 L113
63-900	Large	W35 L127
63-905	X Large	W33 L150
63-906	2X Large	W33 L160
63-907	3X Large	W33 L170
63-908	4X Large	W33 L180



Euro-Med™ Stainless Steel Specula



Vu-Max™ C€

Offered in both Pederson and Graves configurations.

NUMBER	DESCRIPTION	SIZE MM
64-113	Pederson Medium	W22 L114
64-115	Pederson Large	W25 L120
64-218	Graves Medium	W33 L114
64-219	Graves Large	W35 L125
64-225	X Large	W35 L150
64-226	2X Large	W35 L160
64-227	3X Large	W35 L170
64-228	4X Large	W35 L180



Snowman Specula™

Longer blades expose access to the cervix.

NUMBER	DESCRIPTION	SIZE MM
64-107	Graves Blade	W35 L152 C€
64-117	Pederson Blade	W25 L152 C€
64-127	One each/Graves & Pederson	



ExpandaView™ Integrated Speculum and Lateral Retractor

Combines a traditional speculum with a lateral wall retractor.

NUMBER	DESCRIPTION	SIZE MM
64-161	Medium	W22 L109 C€
64-162	Large	W22 L133 C€
64-163	Stainless Steel Instrument Holder	



Auvard ☾

This 2.5-pound weighted speculum for posterior vaginal wall retraction in vaginal procedures.

NUMBER	DESCRIPTION	SIZE MM
64-422	Med (Wide/Short Blade)	W43 L55
64-421	Lg (Narrow/Long Blade)	W35 L90



Weisman Graves Side Opening ☾

Available in two sizes with left or right openings.

NUMBER	DESCRIPTION	SIZE MM
391-667A	Left side open, Med	W33 L102
391-667	Left side open, Lg	W36 L121
391-668A	Right side open, Med	W35 L102
391-668	Right side open, Lg	W36 L120



Pederson Side Opening ☾

The Pederson open-sided specula features either a left or right side opening.

NUMBER	DESCRIPTION	SIZE MM
64-112L	Left side open, Med	W22 L114
64-114L	Left side open, Lg	W25 L120
64-112R	Right side open, Med	W22 L114
64-114R	Right side open, Lg	W25 L120

Euro-Med™ Stainless Steel Specula



Specula Blades and Yoke Dimensions

Use our speculum sizing guide to select the proper tools for your practice. Specula are available in a variety of sizes.

*Approximate size and shape.
See pages 7–8 for a full product selection and model numbers.

Vu-Max™ Yoke – Graves Blades

Medium – W: 33 mm L: 114 mm
Large – W: 35 mm L: 125 mm

Vu-Max Yoke – Pederson Blades

Medium – W: 22 mm L: 114 mm
Large – W: 25 mm L: 120 mm

C. Vu-More™ Yoke – Graves Blades

Mini – W: 34 mm L: 93 mm
Medium – W: 33 mm L: 113 mm
Large – W: 35 mm L: 127 mm

D. Standard Yoke – Pederson Blades

Virgin – W: 13 mm L: 78 mm
Small – W: 13 mm L: 82 mm
Extra Narrow – W: 17 mm L: 108 mm
Medium – W: 22 mm L: 114 mm
Large – W: 25 mm L: 130 mm

Also Available with
Side Opening Yoke

E. Standard Yoke – Graves Blades

Small – W: 19 mm L: 82 mm
Medium – W: 33 mm L: 114 mm
Large – W: 35 mm L: 127 mm

Also Available with
Side Opening Yoke

F. Vu-More Yoke – Graves or Pederson Blade

Snowman™ – W: 35 mm L: 152 mm
Snowman – W: 25 mm L: 152 mm

G. Vu-Max and Vu-More Yokes – Graves Blades

X Large – W: 35 mm L: 150 mm
2X Large – W: 35 mm L: 160 mm
3X Large – W: 35 mm L: 170 mm
4X Large – W: 35 mm L: 180 mm

Euro-Med™ Stainless Steel Lateral Wall Retractors



A.

A. Cer-View™ Lateral Wall Retractor CA CE

NUMBER	DESCRIPTION	SIZE MM
64-310	Medium	W22 L150
64-320	Large	W22 L160






B.

**B. Cer-View™ Lateral
Wall Retractor** €€

NUMBER	DESCRIPTION	SIZE MM
64-307	Extra Large	W32 L190



C.

C. Jackson Vaginal Retractor €€

NUMBER	DESCRIPTION	SIZE MM
64-322	Medium	W40 L98
64-324	Large	W40 L115

Euro-Med™ Stainless Steel Uterine Sounds and Dilators

A. Dil-Os™ Uterine Sound/Dilator

NUMBER	DESCRIPTION	SIZE
64-336	Slight Curve	L12.7"/322 mm
64-337	Strong Curve	L12.5"/317 mm
64-338	Set of two (64-336, 64-337)	



B. Sims Uterine Sound

NUMBER	DESCRIPTION	SIZE
64-600	Sims Uterine Sound Silver-Plated Copper	L12.5"/317
64-601	Sims Uterine Sound Stainless Steel	L12.5"/317

Euro-Med™ Stainless Steel Uterine Sounds and Dilators



Pratt Dilator Set €€

Available individually or in a set of 8.

NUMBER	SIZE
391-830	13 Fr/15 Fr – 41 Fr/43 Fr (8/Set)
391-833	13 Fr/15 Fr
391-834	17 Fr/19 Fr
391-835	21 Fr/23 Fr
391-836	25 Fr/27 Fr
391-837	29 Fr/31 Fr
391-838	33 Fr/35 Fr
391-839	37 Fr/39 Fr
391-840	41 Fr/43 Fr
391-841	45 Fr/47 Fr
391-842	49 Fr/51 Fr
391-843	53 Fr/55 Fr
391-844	57 Fr/59 Fr



Hegar Dilator Set €€

Stainless steel Hegar dilators provide the option of double-ended dilators. Available individually or in a set of 8.

NUMBER	SIZE MM
61207	9 Fr (3 mm)/12 Fr (4 mm)/ 51 Fr (17 mm)/54 Fr (18 mm) (8/Set)
61208	9 Fr (3 mm)/12 Fr (4 mm)
61209	15 Fr (5 mm)/18 Fr (6 mm)
61210	21 Fr (7 mm)/24 Fr (8 mm)
61211	27 Fr (9 mm)/30 Fr (10 mm)
61212	33 Fr (11 mm)/36 Fr (12 mm)
61213	39 Fr (13 mm)/42 Fr (14 mm)
61214	45 Fr (15 mm)/48 Fr (16 mm)
61215	51 Fr (17 mm)/54 Fr (18 mm)



Hank Cervical Dilators €€

Hank cervical dilators available individually or in a set of six ranging from 9 French (3 mm) to 20 French (6.67 mm). The set comes in a custom-fitted and numbered pouch.

NUMBER	SIZE MM
64-370	6 Dilators and Storage Pouch
61201	9 Fr (3 mm)/10 Fr (3.33 mm)
61202	11 Fr (3.67 mm)/12 Fr (4 mm)
61203	13 Fr (4.33 mm)/14 Fr (4.67 mm)
61204	15 Fr (5 mm)/16 Fr (5.33 mm)
61205	17 Fr (5.67 mm)/18 Fr (6 mm)
61206	19 Fr (6.33 mm)/20 Fr (6.67 mm)



Mini Cervical Dilator Set €€

NUMBER	SIZE MM
64-374	5 Dilators, Handle and Storage Pouch



Small Cervical Dilator Set €€

NUMBER	SIZE MM
64-372	6 Dilators and Storage Pouch

Euro-Med™ Stainless Steel Endocervical Curettes



Kevorkian with Basket

NUMBER	LENGTH
66-442	12"/305 mm



Handle Style

Kevorkian without Basket

NUMBER	LENGTH
66-440	12"/305 mm



Townsend™ CE

NUMBER	LENGTH
66-445	12.3"/312 mm



Kevorkian Collector™

NUMBER	LENGTH
66-441	12"/305 mm



London™

NUMBER	LENGTH
64-490	11.7"/297 mm



Handle Style

Endocrette™ CE

NUMBER	LENGTH
66-447	12.3"/312 mm

Euro-Med™ Endometrial Curettes



A. Novak €€

NUMBER	SIZE	LENGTH
64-624	1 mm	9.8"/248 mm
64-626	2 mm	9.8"/248 mm
64-628	3 mm	9.8"/248 mm
64-630	4 mm	9.8"/248 mm

B. Randall €€

NUMBER	SIZE	LENGTH
64-620	3.9 mm	9.9"/252 mm



C. Sims

NUMBER	SIZE	LENGTH
392-110	6 mm	11"/280 mm €€
392-111	7 mm	11"/280 mm €€
392-112	8 mm	11"/280 mm €€
392-113	9 mm	11"/280 mm €€
392-114	11.5 mm	11"/280 mm €€
392-116	17 mm	11"/280 mm

Euro-Med™ Stainless Steel Hooks



Angle Hooks

NUMBER	DESCRIPTION	LENGTH
36-630	Straight Hook	9.8"/249 mm
36-632	Right Hook	9.6"/244 mm



Iris Hooks

NUMBER	LENGTH
36-625	10.3"/262 mm



Marlow Corkscrew Hook CE

NUMBER	LENGTH
67-451	9.3"/236 mm



Double IUD Hook

NUMBER	DESCRIPTION	LENGTH
392-151	10.6"/269 mm	9.8"/249 mm



Emmett Tenaculum Hook CE

NUMBER	LENGTH
66-450	9.8"/249 mm



IUD Hook CE

NUMBER	LENGTH
392-148	10.3"/262 mm

Euro-Med™ Stainless Steel Cervical Tenacula



Braun Tenaculum Forceps €€

NUMBER	LENGTH
64-420	9.9"/252 mm



Schroeder Tenaculum €€

NUMBER	LENGTH
64-410	9.9"/252 mm



Duplay Tenaculum Double Curved €€

NUMBER	LENGTH
61352	11"/280 mm



Goldstein Grasp Cervical Stabilizer €€

NUMBER	LENGTH
64-6026	9.4"/239 mm

Euro-Med™ Stainless Steel Dressing and Tissue Forceps



A. Allis €€

NUMBER	DESCRIPTION	LENGTH
61981	4 x 5 Teeth	6"
61983	5 x 6 Teeth	7.5"
61882	5 x 6 Teeth	10"



B.

B. Dressing Serrated €€

NUMBER	LENGTH
61804	5.6"
61805	6.2"



C.

C. Potts-Smith Tissue €€

NUMBER	DESCRIPTION	LENGTH
61871	1 x 2 Teeth	7"
61872	1 x 2 Teeth	8.2"

Euro-Med™ Stainless Steel Dressing and Tissue Forceps



Potts-Smith Dressing €€

NUMBER	LENGTH
10-202	8.25"



Tissue 1 x 2 €€

NUMBER	DESCRIPTION	LENGTH
61841	1 x 2 Teeth	5.75"
61843	1 x 2 Teeth	7"
61844	1 x 2 Teeth	8"



Adson Dressing Serrated €€

NUMBER	DESCRIPTION	LENGTH
61830	Serrated Tips	4.75"



Adson Tissue 1 x 2 €€

NUMBER	DESCRIPTION	LENGTH
61863	1 x 2 Teeth	4.75"



Adson Tissue 1 x 2 Serrated €€

NUMBER	DESCRIPTION	LENGTH
61864	1 x 2 Serrated Tips	4.75"



Adson-Brown Grasping €€

NUMBER	DESCRIPTION	LENGTH
61866	7 x 7 Teeth	4.75"

Euro-Med™ Stainless Steel Endocervical Specula



Ratchet Handle Speculum €€

Stainless steel endocervical speculum with quick-release lock.

NUMBER	TIP SIZE	LENGTH
66-437	Small 3 mm	10.6"/270 mm
66-435	Medium 4 mm	10.6"/270 mm
66-436	Large 5 mm	10.6"/270 mm



Screw/Knob Handle Speculum €€

The screw/knob handle stainless steel speculum adjusts and maintains position using the threaded thumb screw.

NUMBER	TIP SIZE	LENGTH
66-431	Small 3 mm	10.6"/270 mm
66-430	Medium 4 mm	10.6"/270 mm
66-432	Large 5 mm	10.6"/270 mm



No-Lock Handle €€

This stainless steel endocervical speculum offers single-handed action. The No-Lock Handle allows for adjustment as necessary.

NUMBER	TIP SIZE	LENGTH
66-433	Small 3 mm	10.6"/270 mm
67-435	Medium 4 mm	10.6"/270 mm
66-434	Large 5 mm	10.6"/270 mm

Each Endocervical Speculum Available in Small, Medium or Large Tip

3 mm
Small Size "Thin Line" Tip



4 mm
Medium Size Euro-Med™ Tip



5 mm
Large Size Euro-Med™ Tip



Foerster Sponge Forceps €€

Featuring oval tips with serrations, the Foerster is used to hold sponges during the procedure.

NUMBER	LENGTH
16-182	9.7"

Euro-Med™ Stainless Steel Uterine Forceps



**Bozeman
Uterine Forceps** €€

NUMBER	LENGTH
392-660	10.2"/257 mm



**Kelly
Uterine Forceps** €€

NUMBER	LENGTH
417-131	12.5"/318 mm



**Foerster
Uterine Polyp Forceps** €€

NUMBER	DESCRIPTION	LENGTH
415-140	Straight	10"/254 mm
415-141	Curved	10"/254 mm

Euro-Med™ Stainless Steel Needle Holders



Mayo-Hegar (Wider Jaw) €€

NUMBER	LENGTH
62809	5.5"
62810	6.3"
62815	7.1"
62820	8.1"



Heaney (Curved) €€

NUMBER	LENGTH
62825	8.3"



Crile-Wood (Narrow Jaw) €€

NUMBER	LENGTH
62840	7.1"

Euro-Med™ Stainless Steel Scissors



Operating Scissors (Straight) €€

NUMBER	DESCRIPTION	LENGTH
62002	B/B	5.5"
62008	S/B	5.5"
62014	S/S	5.5"
62004	B/B	6.5"
62005	B/B	7.5"
62010	S/B	6.5"
62011	S/B	7.23"



Spencer Suture Scissors €€

NUMBER	DESCRIPTION	LENGTH
62164	Straight	3.5"



Metzenbaum Scissors (Straight or Curved) €€

NUMBER	DESCRIPTION	LENGTH
62116	Straight	7"
62127	Straight	8"
62120	Straight	9.1"
62122	Straight	11.2"
62131	Curved	5.9"
62118	Curved	7"
62119	Curved	8.1"
62124	Curved	9"
62126	Curved	11"



Mayo Scissors (Straight or Curved) €€

NUMBER	DESCRIPTION	LENGTH
62063	Straight	5.7"
62065	Straight	6.8"
62066	Straight	9"
62067	Curved	5.5"
62069	Curved	6.75"
62070	Curved	9"



Iris Scissors (Straight or Curved) €€

NUMBER	DESCRIPTION	LENGTH
62140	Straight S/S	4.5"
62143	Curved S/S	4.5"



Committed to Our Customers

At CooperSurgical, we strive for the highest level of expertise and knowledge of our products and to develop an ongoing relationship with each of our customers that is built on trust, quality and consistency. We look to provide service that exceeds not only our customers' expectations, but industry standards as well. CooperSurgical's Customer Service Department provides full-service order entry with extended hours to facilitate order placement. As an added convenience we also offer online order entry. Call now to have your Customer Service Representative set up your account today. Once you are registered for online ordering, you will be able to place orders, track shipments, review prior purchase history and schedule orders for future delivery.

CE Available for sale where CE mark is required

CA Available for sale in Canada

CE and CA labeling is subject to change. Contact Customer Service or your local sales representative for the most up-to-date information.

5 Ways to Order

Local Sales Representative

www.coopersurgical.com/pages/salesrepsearch.aspx

Email

U.S.: Orders@coopersurgical.com

Int'l: InternationalCustomerservice@coopersurgical.com

E-Ordering (U.S. Only)

Call 800.243.2974 to set up your account today.

Phone Orders

800.243.2974 or 203.601.5200

(8:00 AM - 8:30 PM EST)

Int'l: +1 203.601.9818

Online

ecommerce.coopersurgical.com



SCAN TO
LEARN MORE

Warranty Information

CooperSurgical, Inc. ("CooperSurgical") warrants that its products offered in this catalog, **other than** single-use products and consumables, will be free from defects in materials and workmanship for one year from the date of purchase.

CooperSurgical warrants that its single-use products and consumables offered in this catalog will, through the earlier of 90 days from the date of purchase of each such product and the expiration date of that product stated on that product's package, be free from defects in materials and workmanship that are discovered prior to or during the single use of that product.

If CooperSurgical determines that a product offered in this catalog fails to conform to the CooperSurgical warranty applicable to that product, CooperSurgical will, as the sole remedy for that failure to conform, repair or replace that product, at CooperSurgical's discretion, free of charge.

To return a product to CooperSurgical, a customer must comply with CooperSurgical's Returned Goods Policy described in this catalog. A customer will not have any remedy if a CooperSurgical product does not conform to the warranty for that product unless that product is returned to CooperSurgical in accordance with that Returned Goods Policy. CooperSurgical will ship returned products that it repaired or replaced under warranty to the customer who returned those products, at CooperSurgical's expense F.O.B. the customer's facility. Under all other circumstances, CooperSurgical will ship returned products to the customer who returned those products at the customer's expense F.O.B. CooperSurgical's facility.

CooperSurgical's warranties, including those in this catalog, do not cover damage caused by misuse, improper care, improper use of chemicals or cleaning methods, loss, theft, servicing or sharpening by non-authorized personnel or negligent or intentional conduct on the part of the owner or user of the product. CooperSurgical's warranties do not cover normal wear and tear or general maintenance, including cleaning and sharpening services required as a result of normal use. Any modifications or changes to a product will void that product's warranty.

THE ABOVE WARRANTIES ARE IN LIEU OF, AND COOPERSURGICAL HEREBY DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, WRITTEN OR ORAL, WITH RESPECT TO COOPERSURGICAL'S PRODUCTS, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO TERMS, CONDITIONS, UNDERSTANDINGS OR AGREEMENTS THAT PURPORT TO MODIFY THE ABOVE WARRANTIES OR THAT MAKE ANY ADDITIONAL WARRANTIES FOR ANY COOPERSURGICAL PRODUCT SHALL HAVE ANY LEGAL EFFECT UNLESS MADE IN WRITING AND SIGNED BY AN AUTHORIZED COOPERSURGICAL CORPORATE OFFICER.

COOPERSURGICAL SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR LOST PROFITS, DAMAGES FROM LOSS OF USE OR LOST DATA, OR INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ITS WARRANTIES OR OTHERWISE FOR ANY CLAIM RELATED TO COOPERSURGICAL'S PRODUCTS, EVEN IF COOPERSURGICAL HAS BEEN ADVISED, KNEW OR SHOULD HAVE KNOWN OF THE POSSIBILITY OF SUCH DAMAGES. COOPERSURGICAL'S LIABILITY WITH RESPECT TO A PRODUCT COVERED BY A WARRANTY OR OTHERWISE SHALL BE LIMITED IN ALL CIRCUMSTANCES TO THE PURCHASE PRICE OF THAT PRODUCT.

Returned Goods Policy

- 1. Goods will be accepted for return for the following reasons:**
 - If shipment was made without Customer's authorization or order
 - If incorrect goods were shipped
 - If defective goods are covered by the standard warranty
- 2. To return goods, you must contact a Customer Service Representative by telephone (1.800.243.2974 or 1.203.601.5200) for a Returned Merchandise Authorization (RMA) number. Items will not be accepted without an RMA number.**

Please have the following information:

 - Reason you wish to return the goods
 - Quantity, description, part number, serial number of the goods
 - Date of Receipt of Order
 - Customer's Purchase Order and the CooperSurgical Invoice Number

All used instruments must be cleaned and sterilized prior to being shipped for return.
- 3. Shipment must be sent prepaid by Customer. Freight collect shipments will not be accepted, and goods will be returned to sender.**
- 4. If Customer intends to return equipment ordered in error, the following restocking charges and terms will apply:**
 - There is no restocking fee if goods are returned within 30 days of shipment
 - There is a restocking fee of 25% if goods are returned before 31-60 days of shipment
 - Goods must be returned unused, in the original carton, and in marketable condition
 - Refurbishing and replacement charges will be added to the restocking charges for damaged or missing items
 - No return after 60 days
 - No refund on sterile, single-use disposable products

**Send All Returns to: CooperSurgical
95 Corporate Drive
Trumbull, CT 06611**



CooperSurgical®

800.243.2974
coopersurgical.com
95 Corporate Drive, Trumbull, CT 06611