



Carter- Thomason CloseSure System®

by  CooperSurgical®



Designed for reproducible, safe, and efficient port site closure

- ✓ Safely and quickly close port site
- ✓ Application across multiple procedure specialties including laparoscopic and robotic procedures (i.e. general, urology, gynecology, and hernia repair)
- ✓ Applicable for a broad range of patients, including obese and bariatric
- ✓ Comprehensive product offering features a broad range of kit configurations and suture passers



Choose the System That's Right for Your Surgical Procedure


The Carter-Thomason Port Closure Portfolio provides a comprehensive offering to ensure consistent suture passings on any body type.

A Carter-Thomason CloseSure System® Suture Passer

Features a ring-style handle for engaging the tip and a recessed hinge jaw needle tip to ensure quick and precise suturing techniques.

B Advanced SureGrasp Suture Passer™

Features an ergonomic plunger mechanism handle, with blunt extendable grasping fingers, that facilitate safe and precise suture capture. Applicable for average, obese, and bariatric populations.



Carter-Thomason CloseSure System®

Number	Description
CTI-512N	(1) 5 mm Pilot Guide 10/12 mm Pilot Guides + Suture Passer CTI-1015N, PG-15S, PG-5/10
CTI-1015N	Carter-Thomason CloseSure System Suture Passer (1) 10/12 mm Pilot Guide (1) 15 mm Pilot Guide, sterile/single-use (5 Systems/Box)
PG-15S	(1) 15 mm sterile/single-use Pilot Guide (5 Sets/Box)
PG-5/10	(1) 5 mm Pilot Guide (1) 10/12 mm Pilot Guide, sterile/single-use (5 Sets/Box)
CTI-SP	(1) Carter-Thomason CloseSure System Suture Passer (5 Passers/Box)

Carter-Thomason CloseSure System XL®

Number	Description
CTXL	Single-use kit includes: (Comes in pack of 5) (1) 10/12 mm Pilot Guide (1) 15 mm Pilot Guide (1) XL CloseSure System Suture Passer

Carter-Thomason II Port Closure System™

Number	Description
CTI-1012P	(1) Suture Passer (1) 10/12 mm Suture Guide (5 Systems/Box)
CTI-15P	(1) Suture Passer (1) 15 mm Suture Guide, (5 Systems/Box)
CTI-1015P	(1) Suture Passer (1) 10/12 mm Suture Guide (1) 15 mm Suture Guide (5 Systems/Box)

Carter-Thomason® SG with SureGrasp Technology

Number	Description
CTI-1015SG	(1) Carter-Thomason SureGrasp Suture Passer (1) 10/12 mm Pilot Guide (1) 15 mm Pilot Guide sterile/single-use (5 Systems/Box)
CTXL-SG	XL Suture Passer (1) 10/12 mm Pilot Guide (1) 15 mm Suture Guide (Box of 3)
CTSG	SureGrasp Suture Passer (1) 10/12 mm Suture Guide (1) 15 mm Suture Guide (Box of 5)