

Anospec[®]

Single-Use Anoscope with Integrated LED Light Source



Disposable, Bright Built-In Light, Batteries Pre-Installed

This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the medical setting.

Anospec is designed to aid in the examination of the anus to screen, diagnose, evaluate, and treat problems of the anal canal. Anospec is ideal for diagnostic and therapeutic anoscopy procedures.

Product Features

- Disposable, single-use device reduces risk of cross-contamination
- Each Anospec comes ready for use with an integrated LED light source and obturator
- No assembly required and no additional parts needed
- Battery-operated LED light source can produce peak light for 20+ minutes
- LED light source generates no damaging heat
- The device's clear material allows for 360 degree visibility



CooperSurgical[®]



Anospec[®]

Ready for Use with an Integrated LED Light Source and Obturator



90 mm x 18 mm Slotted



96 mm x 23 mm Slotted




103 mm x 18 mm Beveled

Ordering Information

NUMBER	DESCRIPTION
C060100	Anospec [®] (90 mm x 18 mm Slotted)
C060110	Anospec [®] (96 mm x 23 mm Slotted)
C060120	Anospec [®] (103 mm x 18 mm Beveled)

Batteries included in each unit. Sold clean, non-sterile

800.243.2974 | 203.601.5200
www.coopersurgical.com



IMPORTANT SAFETY INFORMATION Anospec[®] single-use anoscope is a lighted anal speculum intended for physician use to examine the anal sphincter, anus, rectum, and, using additional devices, to perform various diagnostic and therapeutic procedures.
Warnings and Cautions: Never use with laser equipment, gamma radiation, autoclave sterilization. Used specula should be treated as a biohazardous, infectious material and should be disposed of in a suitable disposable unit in accordance with local regulations. For detailed instructions and risk information, please refer to Instructions for Use. www.coopersurgical.com/anospec-instructions-for-use



CooperSurgical[®]

© 2022 CooperSurgical, Inc. C-US-OBP-000006 Rev. 10/22