



Anospec[®]

Single-Use Anoscope with Integrated LED Light Source



Disposable, Bright Built-In Light, Batteries Pre-Installed

This ready-to-use device eliminates the time and expense of reprocessing and reduces the risk of cross-contamination in the medical setting.

Anospec is designed to aid in the examination of the anus to screen, diagnose, evaluate, and treat problems of the anal canal. Anospec is ideal for diagnostic and therapeutic anoscopy procedures.

Product Features

- Disposable, single-use device reduces risk of cross-contamination
- Each Anospec comes ready for use with an integrated LED light source and obturator
- No assembly required and no additional parts needed
- Battery-operated LED light source can produce peak light for 20+ minutes
- LED light source generates no damaging heat
- The device's clear material allows for 360 degree visibility



CooperSurgical[®]





Anospec[®]

Ready for Use with an Integrated LED Light Source and Obturator



90mm x 18mm Slotted



96mm x 23mm Slotted



103mm x 18mm Beveled

Ordering Information

NUMBER	DESCRIPTION
C060100	Anospec [®] (90mm x 18mm Slotted)
C060110	Anospec [®] (96mm x 23mm Slotted)
C060120	Anospec [®] (103mm x 18mm Beveled)

Batteries included in each unit. Sold clean, non-sterile

800.243.2974 | 203.601.5200
www.coopersurgical.com



CooperSurgical[®]

© 2021 CooperSurgical, Inc.

C-US-OBP-000006

