

Apple-Hunt Secondary Cannulas

- **Uninhibited Abdominal Access**
No bulky external valves, fittings or housings
- **Accepts 3 mm Instruments**
Without reducers
- **Easy Passage of Instruments**
- **No Glare, No Reflection**
Laser reflection is eliminated with specially designed sleeve color and finish
- **Increased Instrument Maneuverability**
Short, low internal and external profile design
- **Two Tip Styles**
Choice of Pyramidal or Conical steel tips
- **Hold Instruments in Place**
Airtight seal clears operative field, holds instruments securely, permitting hands-free exchanges without loss of pneumoperitoneum



Pyramidal Tip




Conical Tip



Airtight Seal

Apple-Hunt Secondary Cannulas



Smooth Cannula/Trocar Transition
Designed for easier sleeve insertion and to lessen the risk of patient injury

Cannula Anchors Securely to Patient
Specially designed, clinically proven fascia-anchoring threads hold cannula securely in place, regardless of number of instrument exchanges

Suturing Through Cannula
Design permits suturing through the cannula

Easy to Work With
Ergonomic design of trocar palm piece makes device easy to handle, convenient to use

5 mm Apple-Hunt Cannula / Trocar			
Number	Description	Number	Description
900-800	5 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar (5/Box)	900-831	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar (3 Pack) (5/Box)
900-801	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar (5/Box)	900-840	5 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar – Extended Length (5/Box)
900-802	5 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with Side Port (5/Box)	900-841	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar – Extended Length (5/Box)
900-803	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar with Side Port (5/Box)	900-842	5 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar – Extended Length (2 pack) (5/Box)
900-820	5 mm Apple-Hunt Secondary Cannula/ Pyramidal Tip Trocar (2 Pack) (5/Box)	900-844	5 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with Side Port – Extended Length (5/Box)
900-821	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar (2 Pack) (5/Box)	900-845	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar with Side Port – Extended Length (5/Box)
900-822	5 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with Side Port (2 Pack) (5/Box)	900-846	5 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with Side Port – Extended Length (2 Pack) (5/Box)
900-823	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar with Side Port (2 Pack) (5/Box)	900-847	5 mm Apple-Hunt Secondary Cannula/Conical Tip Trocar with Side Port – Extended Length (2 Pack) (5/Box)

12 mm Apple-Hunt Cannula / Trocar		10 mm Apple-Hunt Cannula / Trocar		7 mm Apple-Hunt Cannula / Trocar	
Number	Description	Number	Description	Number	Description
900-850	12 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with 5 mm and 10 mm Reducers (5/Box)	900-860	10 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with 5 mm Reducer (5/Box)	900-870	7 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with 5 mm Reducer (5/Box)
900-851	12 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with 5 mm and 10 mm Reducers – Extended Length	900-861	10 mm Apple-Hunt Secondary Cannula/Pyramidal Trocar with 5 mm Reducer – Extended Tip Length (5/Box)	900-871	7/8 mm Apple-Hunt Secondary Cannula/Pyramidal Tip Trocar with 5 mm Reducer – Extended Length (5/Box)
		900-863	10 mm Blunt Tip Trocar/Cannula with 5 mm Reducer (5/Box)		



CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products, which aid clinicians in the management and treatment of commonly seen conditions.

800.243.2974 | 203.601.5200 | coopersurgical.com | 95 Corporate Drive | Trumbull, CT 06611
Also on [YouTube.com/coopersurgical](https://www.youtube.com/coopersurgical)