

Premier Products for Women's Healthcare

Clinic & Practice-Based Products Catalog





We aim to help clinicians strengthen the important relationship between a woman and her OB/Gyn. Our broad portfolio includes products for traditional reproductive medicine, from contraception to fertility and prenatal care, as well as for screening and treatment for urinary incontinence, abnormal uterine bleeding (AUB) and cancer screening.

Clinic and Practice-Based Brands

CooperSurgical® is empowering office-based care with an expanded portfolio of innovative devices and products. Within this catalog, you will find an extensive range of instruments and devices that can be used in your practice to improve patient care.

- Endosee® Advance
- Mara®
- Milex®
- LEEP PRECISION™
- Fischer Cone Biopsy Excisor®
- Pipelle®
- Leisegang™
- ER-Spec®
- Her-Spec®
- Euro-Med™
- H/S Elliptosphere®
- ABBI® Air Bubble Based Infuser
- Paragard® Intrauterine Copper Contraceptive
- Cord Blood RegistrySM
- Wallach®

Table of Contents

| | |
|--|-------|
| Office Hysteroscopy | |
| Endosee® Advance | 6 |
| Water Vapor Ablation System | |
| Mara® | 8-9 |
| Milex® Pessaries | |
| Milex® Pessaries for Prolapse and Incontinence | 10-13 |
| LEEP | |
| LEEP PRECISION™ Integrated System | 14-15 |
| LEEP Electrodes | 16-17 |
| Fischer Cone Biopsy Excisor® | 18 |
| LEEP Procedure Essentials | 19 |
| LEEP System Accessories | 20-21 |
| Euro-Med™ LEEP Instruments | 22-23 |
| Coated Retractors | 24 |
| LEEP Tenacula, Hooks and Rulers | 26-27 |
| Wallach® Electrosurgical System & Accessories | 28-29 |
| Colposcopy | |
| Leisegang OptiK® Series | 30-33 |
| OptiK Model 1 | 30-31 |
| OptiK Model 2 | 32-33 |
| Wallach® Zoom Magnification Colposcopes | 34 |
| Gynecological Specialty Instruments | |
| Lighted Specula and Other Lighted Devices | 36-37 |
| Euro-Med™ Stainless Steel Biopsy Punches | 38-39 |
| Euro-Med™ Rotating Biopsy Punch Handles | 40 |
| Euro-Med™ Standard Biopsy Punches | 41 |
| Euro-Med™ Stainless Steel Specula | 42-45 |
| Stainless Steel Endocervical Specula | 46-47 |
| Stainless Steel Lateral Wall Retractors | 48 |
| Stainless Steel Uterine Sounds and Dilators | 50-51 |
| Stainless Steel Curettes | 52-53 |
| Stainless Steel Hooks | 54 |
| Stainless Steel Cervical Tenacula | 55 |
| Stainless Steel Forceps and Needle Holders | 56-58 |
| Stainless Steel Scissors | 59 |
| Endometrial Samplers | 60-61 |
| Disposable Sound and Dilators | 63 |
| General Examination: Diagnostic Disposables | |
| Dilators | 64 |
| Cervical and Endocervical Samplers | 65 |
| Diagnostic Disposables | 66-69 |
| Cervical Cryosurgery | |
| Leisegang™ LM-900 Cryosurgery System | 70-71 |
| Cryosurgical Systems | 74-75 |
| Fertility and Family Planning | |
| ABBI® | 74 |
| H/S Catheters | 76 |
| Fertility Assessment | 77 |
| Paragard® Intrauterine Copper Contraceptive | 78-79 |
| CBR® Newborn Stem Cell Preservation | 80 |
| Obstetrical Care | |
| Fetal Monitors | 81 |
| Obstetrical Dopplers and Probes | 82-83 |
| Vascular Assessment | |
| Vascular Dopplers and Probes | 84-85 |
| ABI Vascular Systems | 86-89 |
| Dermatologic Cryosurgery | |
| UltraFreeze™ | 91 |
| Index/Customer Service | 92-93 |
| Warranty Information | 93 |

Ⓢ Available for sale where CE mark is required
CA🇺🇸 Available for sale in Canada

Endosee® Advance

In-Office Diagnostic Hysteroscopy Device

Direct visualization system for diagnostic and therapeutic procedures.



Endosee® Advance

| NUMBER | DESCRIPTION | |
|--|---|--------|
| ES9000 | Endosee Advance Display Module, 1/Box | CA☙ CE |
| ESPX5 | Endosee Advance Cannula, Compatible with ES9000 only, 5/Box | CA☙ CE |
| ES-TRAY | Endosee Advance Convenience Kit, 5/Box | |
| ES-TRAY-IV | Endosee Advance IV Convenience Kit, 5/Box | |
| REUSABLE INSTRUMENTS (STAINLESS STEEL) | | |
| ES-BPSY | Biopsy Forceps, 5 Fr, 455 mm, 1/Box | |
| ES-LNGR | Spoon Forceps, Long Serrated, 5 Fr, 455 mm, 1/Box | |
| ES-SCIS | Scissors, Single Action, 5 Fr, 455 mm, 1/Box | |
| SINGLE-USE INSTRUMENTS (STERILE) | | |
| ES-DISP-AGR | Disposable Alligator Grasper Forceps, 5 Fr, 460 mm, 5/Box | |
| ES-DISP-BPY | Disposable Biopsy Cup, 5 Fr, 460 mm, 5/Box | |
| ES-DISP-SCS | Disposable Single Action Scissors, 5 Fr, 460 mm, 5/Box | |



Mara®

Water Vapor Ablation System



ENDOMETRIAL ABLATION DESIGNED TO GENTLY TREAT THE UTERINE LINING, IN THE OFFICE SETTING

- Utilizes the expansive properties of water vapor to conform to each woman's unique uterus, enabling the treatment of more women
- Can significantly reduce heavy menstrual bleeding leading to improved quality of life
- May preserve long-term access to the uterine cavity
- Convenient 2-minute water vapor treatment, as part of a 4-minute procedure*

Risks associated with the use of this procedure are provided in the safety information below.



Mara Water Vapor Ablation System

| NUMBER | DESCRIPTION |
|------------|------------------------------------|
| GEN-16-050 | Mara Water Ablation System Console |
| DDK-16-050 | Mara Water Vapor Probe Kit |


* The water vapor treatment takes 2 minutes, as part of a 4-minute procedure. Total time in office may vary and will include preparation and post-procedure observation.

IMPORTANT SAFETY INFORMATION

Mara Water Vapor Ablation System is indicated to ablate the endometrial lining of the uterus in premenopausal women with menorrhagia due to benign causes for whom childbearing is complete. Pregnancy following the Mara procedure can be dangerous. The Mara procedure is contraindicated for those who have suspected uterine cancer, pre-malignant conditions of the endometrium, endometrial hyperplasia, prior classical cesarean section, transmural myomectomy (including myomectomy performed immediately prior to Mara procedure), taking medications that thin the uterine muscle, uterine length <6cm, prior endometrial ablation, active genital or urinary tract infection, PID, or systemic bacterial infection, IUD currently in place, suspected hydrosalpinx or undiagnosed abnormal vaginal bleeding. Mara is not a sterilization procedure. Rare but serious risks include, but are not limited to, perforation, infection, thermal injury, and serious complications in future pregnancies. Temporary side effects may include cramping, nausea, abdominal pain and distension, vomiting, vaginal infection, and endometritis. Inform patients to contact you if they experience a possible side effect related to Mara. For detailed benefit and risk information, please consult the IFU.

Milex® Pessaries for Prolapse & Incontinence

P = Prolapse I = Incontinence



A. Hodge

For stress urinary incontinence for a woman who has a very small vaginal introitus.

Product # MXPEH (Pessary Only)

I

Hodge

| NUMBER | IN. (A) | MM (A) | IN. (B) | MM (B) |
|---------|---------|--------|---------|--------|
| MXPEH00 | 2 ¾ | 70 | 1 ¾ | 44 |
| MXPEH01 | 3 | 76 | 1 7⁄8 | 48 |
| MXPEH02 | 3 ¼ | 83 | 1 15⁄16 | 49 |
| MXPEH03 | 3 ½ | 89 | 2 | 51 |
| MXPEH04 | 3 5⁄8 | 92 | 2 1⁄8 | 54 |
| MXPEH05 | 3 ¾ | 95 | 2 ¼ | 57 |
| MXPEH06 | 3 7⁄8 | 98 | 2 3⁄8 | 60 |
| MXPEH07 | 4 ¼ | 108 | 2 ½ | 64 |
| MXPEH08 | 4 5⁄8 | 117 | 2 5⁄8 | 67 |
| MXPEH09 | 5 | 127 | 2 ¾ | 70 |

D. Gehrung with Knob/Folding

For cystocele with or without mild rectocele complicated by stress urinary incontinence.

Product # MXPGSK (Pessary Only)

I

E. Gehrung/Folding

Support for cystocele and rectocele. It is also more adaptable in cases of procidentia where the uterus tends to herniate when other pessaries are used.

Product # MXPGS (Pessary Only)

P

F. Cube

Support for 3rd degree prolapse/procidentia. Should be removed nightly.

Product # MXPEC (Pessary Only)

Also available in:
Cube with Drainage Holes
Product # MXPECH (Pessary Only)

P, I

Gehrung

Gehrung w/Knob/Folding

Gehrung Folding

| NUMBER | NUMBER | IN. (A) | MM (A) | IN. (B) | MM (B) | IN. (C) | MM (C) |
|----------|---------|---------|--------|---------|--------|---------|--------|
| MXPGSK00 | MXPGS00 | 1 1⁄16 | 27 | 1 ½ | 38 | 1 3⁄8 | 35 |
| MXPGSK01 | MXPGS01 | 1 1⁄8 | 29 | 1 ¾ | 44 | 1 ½ | 38 |
| MXPGSK02 | MXPGS02 | 1 3⁄8 | 35 | 2 | 51 | 1 5⁄8 | 41 |
| MXPGSK03 | MXPGS03 | 1 ½ | 38 | 2 1⁄8 | 54 | 1 ¾ | 44 |
| MXPGSK04 | MXPGS04 | 1 5⁄8 | 41 | 2 ¼ | 57 | 1 7⁄8 | 48 |
| MXPGSK05 | MXPGS05 | 1 ¾ | 44 | 2 3⁄8 | 60 | 2 | 51 |
| MXPGSK06 | MXPGS06 | 1 7⁄8 | 48 | 2 ½ | 64 | 2 1⁄8 | 54 |
| MXPGSK07 | MXPGS07 | 2 | 51 | 2 7⁄8 | 73 | 2 ¼ | 57 |
| MXPGSK08 | MXPGS08 | 2 1⁄8 | 54 | 3 | 76 | 2 3⁄8 | 60 |
| MXPGSK09 | MXPGS09 | 2 ¼ | 57 | 3 1⁄8 | 79 | 2 ½ | 64 |

Cube

Cube

w/ Holes

| NUMBER | NUMBER | IN. | MM |
|---------|----------|--------|----|
| MXPEC00 | MXPECH00 | 1 | 25 |
| MXPEC01 | MXPECH01 | 1 3⁄16 | 30 |
| MXPEC02 | MXPECH02 | 1 3⁄8 | 35 |
| MXPEC03 | MXPECH03 | 1 ½ | 38 |
| MXPEC04 | MXPECH04 | 1 5⁄8 | 41 |
| MXPEC05 | MXPECH05 | 1 ¾ | 44 |
| MXPEC06 | MXPECH06 | 2 | 51 |
| MXPEC07 | MXPECH07 | 2 ¼ | 57 |

B. Incontinence Dish

For stress urinary incontinence with mild prolapse.

Product # MXPCOND (Pessary Only)

P, I

C. Incontinence Dish w/Support

For stress urinary incontinence with mild prolapse and a mild cystocele.

Product # MXPCONDS (Pessary Only)

P, I

Incontinence Dish

Dish

Dish w/ Support

| NUMBER | NUMBER | IN. | MM |
|-----------|------------|-------|----|
| MXPCOND00 | MXPCONDS00 | 2 1⁄8 | 55 |
| MXPCOND01 | MXPCONDS01 | 2 3⁄8 | 60 |
| MXPCOND02 | MXPCONDS02 | 2 ½ | 65 |
| MXPCOND03 | MXPCONDS03 | 2 ¾ | 70 |
| MXPCOND04 | MXPCONDS04 | 3 | 75 |
| MXPCOND05 | MXPCONDS05 | 3 1⁄8 | 80 |
| MXPCOND06 | MXPCONDS06 | 3 3⁄8 | 85 |

LEEP PRECISION™ Integrated System

Loop Electrosurgical Excision Procedures



Fully integrated and mobile, the LEEP PRECISION System combines three components – Cart, Generator and Smoke Evacuator.

LEEP PRECISION Integrated System

| NUMBER | DESCRIPTION | |
|-----------|--|-----|
| LP-10-120 | LEEP PRECISION Integrated System, 120V | CA* |
| LP-10-220 | LEEP PRECISION Integrated System, 220V International | CE |

LEEP PRECISION Generator

| NUMBER | DESCRIPTION | |
|-----------|--|-----|
| LP-20-120 | LEEP PRECISION Generator, 120V | CA* |
| LP-20-220 | LEEP PRECISION Generator, 220V International | CE |

LEEP PRECISION Smoke Evacuator

| NUMBER | DESCRIPTION | |
|-----------|---|-----|
| LP-30-120 | LEEP PRECISION Smoke Evacuator, 120V, Stand-Alone | CA* |
| LP-30-220 | LEEP PRECISION Smoke Evacuator, 220V, Stand-Alone International | CE |

LEEP PRECISION Cart

| NUMBER | DESCRIPTION | |
|-----------|---------------------|--------|
| LP-40-001 | LEEP PRECISION Cart | CA* CE |

LEEP PRECISION Generator Disposables

| NUMBER | DESCRIPTION | |
|-----------|---|-----|
| LP-50-101 | Hand Switch-Operated Push-Button Cautery Pencils (10/Box) | CA* |
| LP-50-102 | Foot Pedal-Operated Cautery Pencils (10/Box) | CA* |
| LP-50-103 | Foot Pedal Adapter (1/Bag) | CA* |
| LP-50-201 | Dual Dispersive Electrodes with Connector (10/Box) | CA* |
| LP-50-202 | Dual Dispersive Electrodes without Connector (Connect to Reusable Cable) (10/Box) | |
| LP-50-203 | Reusable Cable (1/Bag) | CA* |



LEEP PRECISION Smoke Evacuator Disposables

| NUMBER | DESCRIPTION | |
|--------|--|-----|
| 52560 | Smoke Evacuator Pre-Filter, Single-Patient Use (10/Bag) | CA* |
| 6082 | Smoke Evacuator ULPA Filter (1/Box) | CA* |
| 6083 | Smoke Evacuator Reducer, Single-Patient Use (10/Bag) | CA* |
| 6084 | Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box) | CA* |

LEEP Electrodes

(Shown Actual Size)

Our broad selection of 21 different offerings includes loops, squares, needle and ball tips specifically designed for loop electrosurgical excision procedures.



Large-Radius Loop Electrodes ≥ 2 cm Wide

| NUMBER | DESCRIPTION |
|--------|---|
| R2510 | 2.5 cm W x 1.0 cm D • 12 cm shaft (5/Box) |
| R2020 | 2.0 cm W x 2.0 cm D • 12 cm shaft (5/Box) |
| R2015 | 2.0 cm W x 1.5 cm D • 12 cm shaft (5/Box) |
| R2010 | 2.0 cm W x 1.0 cm D • 12 cm shaft (5/Box) |
| R2008 | 2.0 cm W x 0.8 cm D • 12 cm shaft (5/Box) |

Medium-Radius Loop Electrodes <2 cm Wide

| NUMBER | DESCRIPTION |
|--------|---|
| R1515 | 1.5 cm W x 1.5 cm D • 12 cm shaft (5/Box) |
| R1507 | 1.5 cm W x 0.7 cm D • 12 cm shaft (5/Box) |
| R1505 | 1.5 cm W x 0.5 cm D • 12 cm shaft (5/Box) |
| R1010 | 1.0 cm W x 1.0 cm D • 12 cm shaft (5/Box) |
| R1007 | 1.0 cm W x 0.7 cm D • 12 cm shaft (5/Box) |

Square Electrodes

| NUMBER | DESCRIPTION |
|--------|---|
| S1010 | 1.0 cm W x 1.0 cm D 12 cm shaft (5/Box) |
| S1007 | 1.0 cm W x 0.7 cm D 5.5 cm shaft 15° bend (5/Box) |
| S1004 | 1.0 cm W x 0.4 cm D 5.5 cm shaft 15° bend (5/Box) |
| V1004 | 1.0 cm W x 0.4 cm D 8 cm shaft (5/Box) |

Needle Electrodes

| NUMBER | DESCRIPTION |
|--------|--|
| N1200 | 0.8 mm D x 16 mm 12 cm shaft (5/Box) |
| N0550 | 0.8 mm D x 16 mm 5.5 cm shaft (5/Box) |
| N2512 | 0.25 mm D x 8 mm 12 cm shaft (5/Box) |
| N2555 | 0.25 mm D x 5 mm 5.5 cm shaft (5/Box) |

Ball Electrodes

| NUMBER | DESCRIPTION |
|--------|--------------------------------|
| B0512 | 5 mm D 12 cm shaft (5/Box) |
| B0555 | 5 mm D 5.5 cm shaft (5/Box) |
| B0312 | 3 mm D 12 cm shaft (5/Box) |


Fischer Cone Biopsy Excisor®

For Cone Biopsy Procedures

The Fischer Cone Biopsy Excisor (FCBE) comes in a variety of sizes and shapes.

Fischer Cone Biopsy Excisor CA★ CE

| NUMBER | DESCRIPTION |
|---------|--|
| 900-150 | Fischer Cone Biopsy Excisor, Small (5/Box) |
| 900-151 | Fischer Cone Biopsy Excisor, Medium (5/Box) |
| 900-152 | Fischer Cone Biopsy Excisor, Large (5/Box) |
| 900-154 | Fischer Cone Biopsy Excisor, Medium Extended (5/Box) |
| 900-155 | Fischer Cone Biopsy Excisor, Large Shallow (5/Box) |
| 900-156 | Fischer Cone Biopsy Excisor, Mixed Set (S, M, ME, L & LS) 1 each (5/Box) |
| 900-157 | Fischer Cone Biopsy Excisor, Small Wide-Angle (5/Box) |
| 900-158 | Fischer Cone Biopsy Excisor, Small Wide-Angle Extended (5/Box) |



LEEP Procedure Essentials

Potocky Needle®

27-gauge Potocky Needle. The Potocky Needle attaches to a “dental style” cartridge syringe, while the Potocky Needle with Luer Lock is used with a standard medical syringe.

| NUMBER | DESCRIPTION |
|--------|---|
| 6066 | Potocky Needles (6/Box) |
| 6066-4 | Potocky Needles – 4 Boxes of 6 (24/Case) |
| 6166 | Potocky Needles with Luer Lock (6/Box) |
| 6166-4 | Potocky Needles with Luer Lock – 4 Boxes of 6 (24/Case) |



Needle Extension

| NUMBER | DESCRIPTION |
|---------|----------------------------------|
| 433-132 | Needle Extension • 10 cm (1/Box) |
| 433-133 | Needle Extension • 13 cm (1/Box) |
| 433-134 | Needle Extension • 15 cm (1/Box) |



LEEP System Accessories



Foot Pedal, LEEP System 1000

Featuring pneumatic, air-operated control, the foot pedal does not use electric current.

| NUMBER | DESCRIPTION |
|--------|---------------------------------|
| 6032 | Foot Pedal for LEEP System 1000 |



Handpiece, Single-Patient Use CA*

The sterile, single-patient use handpiece works with any of our disposable electrodes. The wand is operated with a footswitch.

| NUMBER | DESCRIPTION |
|--------|--|
| 6040 | Handpiece, Single-Patient Use (10/Box) |

Patient Return Pad Adapter LEEP System 1000



This reusable replacement adapter connects the single-use Patient Return Pad to the LEEP System 1000.


| NUMBER | DESCRIPTION |
|--------|--|
| 23635 | Patient Return Pad Adapter for LEEP System 1000 (for use with 6050P1 pads) |



Return Pad, Single-Patient Use CA*

Non-sterile patient return pads.

| NUMBER | DESCRIPTION |
|--------|--|
| 6050P1 | LEEP Return Pad, Single-Patient Use (10/Box) |



Disposable Smoke Evacuator (DSE) Tubing

Non-sterile single-use, DSE tubing is designed to work specifically with the LEEP rose-coated specula. The 9-inch DSE tubing is flexible and noncombustible and contains no harmful PVC.

Also available as a set with 6 feet of tubing. See Tubing Set, Single-Patient Use, product item number 6084.

| NUMBER | DESCRIPTION |
|--------|---|
| 6150 | Disposable Smoke Evacuator Tubing – 9" straw (50/Bag) |
| 34371 | LEEP System 1000 Manual (No charge) |

LEEP Accessories



Tubing Set, Single-Patient Use CA*

Each set includes 6 feet of sterile disposable tubing and a DSE (Disposable Smoke Evacuator) tube. This set completes the connection between the reducer and the specula.

| NUMBER | DESCRIPTION |
|--------|--|
| 6084 | Smoke Evacuator Tubing – 6 ft. long, 1/2 in. OD (Sterile, Single-Patient Use) (10/Box) |



Reducer CA*

Disposable, single-patient use reducer connects the Pre-Filter (52560) with the Patient Tubing Set (6084) on the Smoke Evacuator.

| NUMBER | DESCRIPTION |
|--------|----------------------------------|
| 6083 | Smoke Evacuator Reducer (10/Box) |



Pre-Filter CA*

Replacement filters for the CooperSurgical Smoke Evacuator. Used to eliminate smoke odor.

| NUMBER | DESCRIPTION |
|--------|-------------------------------------|
| 52560 | Smoke Evacuator Pre-Filter (10/Box) |



ULPA Filter CA*

Used with the Smoke Evacuator.

| NUMBER | DESCRIPTION |
|--------|-------------------------------------|
| 6082 | Smoke Evacuator ULPA Filter (1/Box) |

Compatible Disposable Accessories

See your local sales representative for more information.

| NUMBER | DESCRIPTION |
|---------|--|
| 395-210 | Electrode Tip Cleanser, Disposable (24/Box) |
| 395-302 | Conductive Electrode Pads, Disposable for ESU (25/Box) |
| 395-305 | Conductive Electrode Pads, Dual Dispersive, Disposable (25/Box) Plus II Gen. |
| 395-310 | Hand-Switch Pencil and Cord, Disposable (25/Box) |
| 72308 | Dual Section Grounding Pads (10/Box) |
| 72102 | Filter/Hose Set for Smoke Evacuator (Pack of 3) |
| 72311 | Patient Return Cable |
| 22192 | LEEP 6000 Return Pad |
| 22191 | Foot Pedal |

Euro-Med™ LEEP Instruments

These nonconductive, fully autoclavable LEEP instruments are offered in a variety of styles to facilitate cervical tissue excision and removal.



Pederson LEEP Specula

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| F223 | Small | W13 L81 |
| F221 | Medium | W22 L114 |
| F220 | Large | W25 L120 |



Vu-More™ LEEP Specula

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| F217 | Medium | W33 L113 |
| F215 | Large | W35 L127 |



Vu-Max™ LEEP Specula

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-----------------|----------|
| F222 | Pederson Medium | W22 L114 |
| F224 | Pederson Large | W25 L120 |
| F218 | Graves Medium | W33 L113 |
| F219 | Graves Large | W35 L127 |



Graves LEEP Specula

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| F205 | Medium | W33 L113 |
| F200 | Large | W35 L127 |



Snowman Specula™

| NUMBER | DESCRIPTION | SIZE MM |
|--------|------------------------------|----------|
| F107 | Graves Blade | W35 L152 |
| F117 | Pederson Blade | W25 L152 |
| F127 | One each/Graves and Pederson | |



Mini Vu-More™ LEEP Specula

| NUMBER | SIZE MM |
|--------|---------|
| F250 | W34 L93 |



Open-Sided LEEP Specula

| NUMBER | SIZE MM |
|--------|----------|
| F210 | W35 L107 |



Auvard Surgical Specula

| NUMBER | DESCRIPTION | SIZE MM |
|--------|----------------------------|---------|
| F230 | Medium (Wide/ Short Blade) | W43 L55 |
| F232 | Large (Narrow/Long Blade) | W35 L90 |




Pederson LEEP Specula (with Integrated Smoke Tube)

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| 61521 | Small | W16 L84 |
| 61518 | Medium | W22 L114 |
| 61519 | Large | W26 L121 |




Graves LEEP Specula (with Integrated Smoke Tube)

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| 61513 | Small | W22 L87 |
| 61514 | Medium | W33 L113 |
| 61515 | Large | W35 L125 |



Euro-Med™ LEEP Instruments



Tru-View™ LEEP
Lateral Wall Retractor

| NUMBER | DESCRIPTION | SIZE MM |
|--------|------------------------|----------|
| F317 | Lateral Wall Retractor | W32 L190 |



Jackson LEEP Retractor

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| F235 | Medium | W40 L89 |
| F234 | Large | W40 L108 |



Cer-View™ LEEP
Lateral Wall Retractor

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| F410 | Medium | W22 L142 |
| F420 | Large | W22 L155 |




Euro-Med™ LEEP Tenacula, Hooks and Rulers



LEEP Ring Forceps

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F920 | 220 mm |



LEEP Tissue Forceps

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F700 | 210 mm |



LEEP Campion Forceps

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F925 | 220 mm |



LEEP Single-Tooth Tenaculum

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F900 | 240 mm |



LEEP Emmett Tenaculum

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F805 | 250 mm |



LEEP Two-Prong Hook

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F802 | 260 mm |



LEEP Three-Prong Hook

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F803 | 260 mm |



LEEP Iris Hook

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F800 | 263 mm |



LEEP Hook

| NUMBER | INSTRUMENT LENGTH | |
|--------|-------------------|--------|
| F810 | Straight | 250 mm |
| F812 | Left-Angled | 250 mm |
| F813 | Right-Angled | 250 mm |



LEEP Depth Ruler

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F600 | 240 mm |



LEEP Width Ruler

| NUMBER | INSTRUMENT LENGTH |
|--------|-------------------|
| F610 | 210 mm |

LEEP Endocervical Specula



Large


| NUMBER | TIP SIZE | LENGTH |
|--------|----------|--------|
| F500 | 5 mm | 273 mm |



Medium


| NUMBER | TIP SIZE | LENGTH |
|--------|----------|--------|
| F510 | 4 mm | 273 mm |

Electrosurgical System Accessories for Wallach® Quantum 2000™



Foot Pedal–Operated Pencils
Sterile, single-use only. For use with Quantum 2000™ and other comparable electrosurgical generators. Reusable adapter required. Use with any standard disposable electrodes.

| NUMBER | DESCRIPTION |
|--------|---|
| 909028 | Foot Pedal–Operated Pencils, Sterile (25 per Box) |



Disposable Split Grounding Pads (with Connector)
For use with Quantum 2000™.

| NUMBER | DESCRIPTION |
|--------|---|
| 909078 | Disposable Split Grounding Pads with Connector (10 per Box) |



EndoBlock Needle (for Intra-Cervical Block)
Sterile, single-use only, reinforced with an extended shaft. Threaded hub needle for dental style cartridge syringe. Luer lock for use with medical syringe.

| NUMBER | DESCRIPTION |
|----------|---|
| 920017 | EndoBlock Needle Luer Lock (6 per Box) |
| 920017-1 | EndoBlock Needle Luer Lock (60 per Box) |
| 920021 | EndoBlock Needle Threaded Hub (6 per Box) |



Disposable Smoke Evacuator Tubing with Attached Speculum Tubing
Sterile, single-use six-foot-long tubing connects to the smoke evacuator filter assembly.

| NUMBER | DESCRIPTION |
|--------|---|
| 920002 | Disposable Smoke Evacuator Tubing with Speculum Tubing (10 per Box) |



Endocervical Block Kit (for Intra-Cervical Block)
Single-use only. Reinforced long 18 gauge EndoBlock needle for drawing anesthetic. Disposable syringe with finger grip.

| NUMBER | DESCRIPTION |
|--------|--|
| 920009 | Endocervical Block Kit (10 Kits per Box) |



Disposable Smoke Evacuator Speculum Tubing (Non-Sterile)

| NUMBER | DESCRIPTION |
|--------|---|
| 920001 | Disposable Speculum Tubing (50 per Box) |



Leisegang OptiK® Series


Optical Devices Designed for Non Contact Visualization of the Tissues of the Cervix, Vulva, Vagina, And Anogenital Areas




OptiK Model 1



Four Base Choices:



Tilt Base
OPTIK1-00



Upright Rolling Base
OPTIK1-01



Swing-Arm Rolling Base
OPTIK1-02



Balance-O-Matic Rolling Base
OPTIK1-03

LEISEGANG OptiK® Model 1

| NUMBER | DESCRIPTION |
|-----------|--|
| OPTIK1-00 | Leisegang OptiK Model 1 Tilt Colposcope (1DL LED): 3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, Tilt Base |
| OPTIK1-01 | Same as OPTIK1-00 except with Upright Rolling Base |
| OPTIK1-02 | Same as OPTIK1-00 except with medium-weight Swing-Arm Rolling Base |
| OPTIK1-03 | Same as OPTIK1-00 except with heavy-weight Balance-O-Matic Rolling Base |

OptiK Model 2


OptiK Model 2 with a Digital SLR Camera



OptiK Model 2 with an Integrated Camera

Three Base Choices:

OptiK Model 2 is available in three base configurations including Upright, Swing-Arm and Balance-O-Matic. The OptiK Model 2 features a working distance of 300 mm, field of view of 19 mm to 76 mm, and depth of field of 2 mm to 8 mm.



Upright Rolling Base
OPTIK2-03-002



Swing-Arm Rolling Base
OPTIK2-04-002



Balance-O-Matic Rolling Base
OPTIK2-05-002

Leisegang OptiK® Model 2 Sales Configuration Options:

LEISEGANG OptiK® Model 2, Base and dSLR Camera CA €

| NUMBER | DESCRIPTION |
|---------------|---|
| OPTIK2-03-002 | 3-Step magnification, 300 mm working distance, LED light, 24,000 Lux, 1.0" CCD optimized, medium-weight Upright Rolling Base, T2 adapter/tubus, Canon camera kit, battery, charger, USB wire, A/C adapter, wired remote and memory card |
| OPTIK2-04-002 | Same as OPTIK2-03-002 except with Swing Arm Rolling Base |
| OPTIK2-05-002 | Same as OPTIK2-03-002 except with heavy-weight Balance-O-Matic Rolling Base |


OptiK Model 2 Accessories

Cameras and Monitors

| NUMBER | DESCRIPTION |
|---------|---|
| 53576 | Canon® Rebel Camera Kit for OptiK Model 2 1.0" CCD only |
| 23676FG | Foot Pedal for Canon Rebel Camera |

Wallach® Zoom Magnification Colposcopes

Continuously Variable Nikon® 4.5 – 20X Zoom,
Trulight LED Lighting



ZoomScope®
with Trulight™ CA★

| NUMBER | DESCRIPTION |
|---------------|--|
| 906043-SP-3 | ZoomScope Overhead arm colposcope Quantum series. Smart purchase package: includes a set of sterilizable handles, a set of sterilizable knob covers for zoom and focus, dust cover; 3-leg base and seven year warranty |
| 906043-SP-4 | ZoomScope with Trulight and 4-Leg Base |
| 906043-40TU-5 | ZoomScope Trulight Video, Digital USB Camera |
| 906057T | ZoomStar with Trulight |
| 906057-40TU | ZoomStar with Trulight Video |



Lighted Specula and Other Lighted Devices



ER-Spec®

Single-use vaginal speculum with integrated LED light source.

| NUMBER | DESCRIPTION |
|-----------|-----------------------------|
| C020130-1 | ER-Spec® (Extra-Small) |
| C020100-1 | ER-Spec® (Small) |
| C020110-1 | ER-Spec® (Medium) |
| C020120-1 | ER-Spec® (Large) |
| C020140-1 | ER-Spec® (Extra-Large) |
| C020300 | ER-Spec® (Small - Sterile) |
| C020310 | ER-Spec® (Medium - Sterile) |
| C020320 | ER-Spec® (Large - Sterile) |



Her-Spec®

Single-use click-free vaginal speculum with integrated LED light source.

| NUMBER | DESCRIPTION |
|---------|---|
| C020500 | Her-Spec® Small, Non-Sterile (1 Box of 16) |
| C020510 | Her-Spec® Medium, Non-Sterile (1 Box of 16) |
| C020520 | Her-Spec® Large, Non-Sterile (1 Box of 14) |
| C020600 | Her-Spec® Small, Sterile (1 Box of 10) |
| C020610 | Her-Spec® Medium, Sterile (1 Box of 8) |
| C020620 | Her-Spec® Large, Sterile (1 Box of 8) |

Accessories

| NUMBER | DESCRIPTION |
|---------|---|
| C040100 | OfficePack® for Hysteroscopy (1 Box of 5) |
| C050100 | Fluid Management Pack (1 Box of 5) |
| C070110 | IUD Pack (1 Box of 10) |



OfficeSpec®

Single-use side-opening vaginal speculum with integrated LED light source.

| NUMBER | DESCRIPTION |
|---------|-----------------------------------|
| C020001 | OfficeSpec® (Side-Opening Medium) |



Anospec®

Designed to aid in the examination of the anus to screen, diagnose, evaluate, and treat problems of the anal canal.

| NUMBER | DESCRIPTION |
|---------|---|
| C060100 | Anospec®, 90 mm x 18 mm slotted, Non-Sterile (1 Box of 20) |
| C060110 | Anospec®, 96 mm x 23 mm slotted, Non-Sterile (1 Box of 16) |
| C060120 | Anospec®, 103 mm x 18 mm beveled, Non-Sterile (1 Box of 20) |

Euro-Med™ Stainless Steel Biopsy Punches

Blue markings indicate bite shape.

A. Tischler-Morgan™ Biopsy Punch

The standard bite incorporates a single tooth on the upper and lower jaws that draws tissue in for sample retention.

| NUMBER | LENGTH | BITE |
|---------|--------|-------------|
| 64-689 | 21 cm | 7 mm x 3 mm |
| 64-689L | 25 cm | 7 mm x 3 mm |

B. Baby Tischler Biopsy Punch

Small bite incorporates a tooth permitting sampling at the endocervical canal.

| NUMBER | LENGTH | BITE |
|---------|--------|-----------------|
| 64-691 | 21 cm | 4.1 mm x 2.3 mm |
| 64-691L | 26 cm | 4.1 mm x 2.3 mm |

C. Tischler-Kevorkian™ Biopsy Punch

Large rectangular bite includes upper and lower teeth.

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-679 | 22 cm | 9 mm x 3 mm |

D. Kevorkian-Pacific™ Biopsy Punch

A slightly smaller bite and additional anchoring teeth on the bottom jaw.

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-675 | 22 cm | 9 mm x 3 mm |

E. Mini-Townsend™ Biopsy Punch

| NUMBER | LENGTH | BITE |
|--------|--------|---------------|
| 64-690 | 21 cm | 4 mm x 2.3 mm |

F. Coppleson Biopsy Punch



The unique jaw shape yields a distinct sample orientation for histologic analysis.

| NUMBER | LENGTH | BITE |
|--------|--------|---------------|
| 64-660 | 22 cm | 7.5 mm x 3 mm |

G. Baggish Biopsy Punch

This punch provides an oval-shaped bite.
The jaws have two bottom teeth and a serrated bed.

| NUMBER | LENGTH | BITE |
|--------|--------|-----------------|
| 64-695 | 22 cm | 6.5 mm x 3.5 mm |



H. Burke™ Biopsy Punch

The design features an open handle and an oval-shaped bite for small samples.

| NUMBER | LENGTH | BITE |
|--------|--------|---------------|
| 64-677 | 22 cm | 5.5 mm x 3 mm |

I. Euro-Med™ Keyes Vulvar Punch Set and Singles


This set includes six stainless steel interchangeable tips ranging from 3 mm to 8 mm in diameter, and a universal handle. The individual tips can be placed on the threaded handle as required.

| NUMBER | LENGTH |
|--------|--|
| 68-180 | Set includes: five tips, storage tray and wooden box |

Euro-Med™ Rotating Biopsy Punch Handles

A 360-Degree Rotating Shaft for Sample Collection

These rotating punch tips fit either the Euro-Med™ closed rotating handle or the Townsend open rotating handle.



Tischler-Morgan™ Tip €€

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-484 | 23 cm | 7 mm x 3 mm |

Mini Tip Down €€

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-656 | 23 cm | 7 mm x 3 mm |

Tischler Tip Up €€

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-751 | 23 cm | 7 mm x 3 mm |

Mini-Townsend™ €€

| NUMBER | LENGTH | BITE |
|--------|--------|---------------|
| 64-460 | 23 cm | 4 mm x 2.3 mm |

Tischler Tip Down €€

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-756 | 23 cm | 7 mm x 3 mm |

Oval Tip €€

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-489 | 23 cm | 5 mm x 3 mm |

Mini Tip Up €€

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-651 | 23 cm | 7 mm x 3 mm |

Kevorkian Tip €€

| NUMBER | LENGTH | BITE |
|--------|--------|-------------|
| 64-488 | 23 cm | 8 mm x 3 mm |

Euro-Med™ Standard Biopsy Punches



Eppendorfer Tip €€

This oval cutting jaw obtains an 8 mm x 4.5 mm cervical specimen. A basket in the jaw helps retain the biopsy sample.

| NUMBER | LENGTH | BITE |
|--------|--------|---------------|
| 64-687 | 21 cm | 8 mm x 4.5 mm |

Wittner Tip €€

The angle of the jaw enables extra access to areas of the cervix. The large teeth and strong bite provide a maximum tissue sample area.

| NUMBER | LENGTH | BITE |
|--------|---------|-------------|
| 64-681 | 21.5 cm | 8 mm x 4 mm |

Kevorkian Tip €€

This punch facilitates simple insertion and sampling of the cervical canal.

| NUMBER | LENGTH | BITE |
|--------|---------|-------------|
| 64-685 | 26.5 cm | 9 mm x 3 mm |

Schubert Tip €€

Features a large bite and an angled handle.


| NUMBER | LENGTH | BITE |
|--------|---------|---------------|
| 64-670 | 21.5 cm | 8 mm x 4.5 mm |

Euro-Med™ Stainless Steel Specula

A. Graves C€

Available with standard yoke and an open-sided configuration. The Graves blade is also available in longer lengths with a Vu-More™ or Vu-Max™ yoke.

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| 64-100 | Small | W19 L82 |
| 64-102 | Medium | W33 L114 |
| 64-104 | Large | W35 L127 |
| 64-108 | Open-sided | W35 L100 |



A.

B. Pederson

Pederson vaginal speculum is a double-bladed instrument with a screw mechanism to keep the blades open once inserted into the vaginal canal. The speculum is designed to aid in the visualization of the vaginal canal for diagnostic and/or procedural purposes.

| NUMBER | DESCRIPTION | SIZE MM |
|---------|--------------|-------------|
| 391-671 | Virgin | W13 L78 |
| 64-110 | Small | W13 L82 C€ |
| 391-674 | Extra-narrow | W17 L108 C€ |
| 64-112 | Medium | W22 L114 C€ |
| 64-114 | Large | W25 L130 C€ |




B.

C. Vu-More™ C€


The Vu-More provides a 25 percent larger viewing area at the yoke.

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| 64-130 | Mini | W34 L93 |
| 63-902 | Medium | W33 L113 |
| 63-900 | Large | W35 L127 |
| 63-905 | X Large | W33 L150 |
| 63-906 | 2X Large | W33 L160 |
| 63-907 | 3X Large | W33 L170 |
| 63-908 | 4X Large | W33 L180 |

New Longer Length



C.



D.

D. Vu-Max™ C€

Offered in both Pederson and Graves configurations.

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-----------------|----------|
| 64-113 | Pederson Medium | W22 L114 |
| 64-115 | Pederson Large | W25 L120 |
| 64-218 | Graves Medium | W33 L114 |
| 64-219 | Graves Large | W35 L125 |
| 64-225 | X Large | W35 L150 |
| 64-226 | 2X Large | W35 L160 |
| 64-227 | 3X Large | W35 L170 |
| 64-228 | 4X Large | W35 L180 |

New Longer Length



E.

E. Snowman Specula™

Longer blades expose access to the cervix.

| NUMBER | DESCRIPTION | SIZE MM |
|--------|----------------------------|-------------|
| 64-107 | Graves Blade | W35 L152 C€ |
| 64-117 | Pederson Blade | W25 L152 C€ |
| 64-127 | One each/Graves & Pederson | |



F.

F. ExpandaView™ Integrated Speculum and Lateral Retractor

Combines a traditional speculum with a lateral wall retractor.

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-----------------------------------|-------------|
| 64-161 | Medium | W22 L109 C€ |
| 64-162 | Large | W22 L133 C€ |
| 64-163 | Stainless Steel Instrument Holder | |



G.

G. Auvard C€

This 2.5-pound weighted speculum for posterior vaginal wall retraction in vaginal procedures.

| NUMBER | DESCRIPTION | SIZE MM |
|--------|------------------------|---------|
| 64-422 | Med (Wide/Short Blade) | W43 L55 |
| 64-421 | Lg (Narrow/Long Blade) | W35 L90 |



H.

H. Weisman Graves Side Opening C€

Available in two sizes with left or right openings.

| NUMBER | DESCRIPTION | SIZE MM |
|----------|----------------------|----------|
| 391-667A | Left side open, Med | W33 L102 |
| 391-667 | Left side open, Lg | W36 L121 |
| 391-668A | Right side open, Med | W35 L102 |
| 391-668 | Right side open, Lg | W36 L120 |



I.

I. Pederson Side Opening C€

The Pederson open-sided specula features either a left or right side opening.

| NUMBER | DESCRIPTION | SIZE MM |
|---------|----------------------|----------|
| 64-112L | Left side open, Med | W22 L114 |
| 64-114L | Left side open, Lg | W25 L120 |
| 64-112R | Right side open, Med | W22 L114 |
| 64-114R | Right side open, Lg | W25 L120 |

(W=Distal End) (L=Penetration Length)

Euro-Med™ Stainless Steel Specula



Specula Blades and Yoke Dimensions

Use our speculum sizing guide to select the proper tools for your practice. Specula are available in a variety of sizes.

See pages 44-45 for a full product selection and model numbers.

Vu-Max™ Yoke – Graves Blades

Medium – W: 33 mm L: 114 mm
Large – W: 35 mm L: 125 mm

Vu-Max Yoke – Pederson Blades

Medium – W: 22 mm L: 114 mm
Large – W: 25 mm L: 120 mm

C. Vu-More™ Yoke – Graves Blades

Mini – W: 34 mm L: 93 mm
Medium – W: 33 mm L: 113 mm
Large – W: 35 mm L: 127 mm

D. Standard Yoke – Pederson Blades

Virgin – W: 13 mm L: 78 mm
Small – W: 13 mm L: 82 mm
Extra Narrow – W: 17 mm L: 108 mm
Medium – W: 22 mm L: 114 mm
Large – W: 25 mm L: 130 mm

Also Available with
Side Opening Yoke

E. Standard Yoke – Graves Blades

Small – W: 19 mm L: 82 mm
Medium – W: 33 mm L: 114 mm
Large – W: 35 mm L: 127 mm

Also Available with
Side Opening Yoke

F. Vu-More Yoke – Graves or Pederson Blade

Snowman™ – W: 35 mm L: 152 mm
Snowman – W: 25 mm L: 152 mm

G. Vu-Max and Vu-More Yokes – Graves Blades

X Large – W: 35 mm L: 150 mm
2X Large – W: 35 mm L: 160 mm
3X Large – W: 35 mm L: 170 mm
4X Large – W: 35 mm L: 180 mm

Euro-Med™ Stainless Steel Endocervical Specula



Ratchet Handle Speculum ☾

Stainless steel endocervical speculum with quick-release lock.

| NUMBER | TIP SIZE | LENGTH |
|--------|-------------|--------------|
| 66-437 | Small 3 mm | 10.6"/270 mm |
| 66-435 | Medium 4 mm | 10.6"/270 mm |
| 66-436 | Large 5 mm | 10.6"/270 mm |

Screw/Knob Handle Speculum ☾

The screw/knob handle stainless steel speculum adjusts and maintains position using the threaded thumb screw.


| NUMBER | TIP SIZE | LENGTH |
|--------|-------------|--------------|
| 66-431 | Small 3 mm | 10.6"/270 mm |
| 66-430 | Medium 4 mm | 10.6"/270 mm |
| 66-432 | Large 5 mm | 10.6"/270 mm |

No-Lock Handle ☾

This stainless steel endocervical speculum offers singlehanded action. The No-Lock Handle allows for adjustment as necessary.

| NUMBER | TIP SIZE | LENGTH |
|--------|-------------|--------------|
| 66-433 | Small 3 mm | 10.6"/270 mm |
| 67-435 | Medium 4 mm | 10.6"/270 mm |
| 66-434 | Large 5 mm | 10.6"/270 mm |

Each Endocervical Speculum Available in Small, Medium or Large Tip



Foerster Sponge Forceps ☾

Featuring oval tips with serrations, the Foerster is used to hold sponges during the procedure.

| NUMBER | LENGTH |
|--------|--------|
| 16-182 | 9.7" |

Euro-Med™ Stainless Steel Lateral Wall Retractors



Cer-View™ Lateral Wall Retractor CA* CE

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| 64-310 | Medium | W22 L150 |
| 64-320 | Large | W22 L160 |

Tru-View™ Lateral Wall Retractor CE

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| 64-307 | Extra Large | W32 L190 |

Jackson Vaginal Retractor CE

| NUMBER | DESCRIPTION | SIZE MM |
|--------|-------------|----------|
| 64-322 | Medium | W40 L98 |
| 64-324 | Large | W40 L115 |



Euro-Med™ Stainless Steel Uterine Sounds and Dilators

Dil-Os™ Uterine Sound/Dilator

| NUMBER | DESCRIPTION | SIZE | |
|--------|-----------------------------|---------------|----|
| 64-336 | Slight Curve | L12.7"/322 mm | €€ |
| 64-337 | Strong Curve | L12.5"/317 mm | €€ |
| 64-338 | Set of two (64-336, 64-337) | | |

Sims Uterine Sound €€

| NUMBER | DESCRIPTION | SIZE |
|--------|--|---------------|
| 64-601 | Sims Uterine Sound – Stainless Steel | L12.5"/317 mm |
| 64-600 | Sims Uterine Sound – Silver-Plated Copper* | L12.5"/317 mm |

* Silver-plated copper may tarnish.



Pratt Dilator Set €€

Available individually or in a set of 8.


| NUMBER | SIZE | NUMBER | SIZE |
|---------|-----------------------------------|---------|-------------|
| 391-830 | 13 Fr/15 Fr – 41 Fr/43 Fr (8/Set) | 391-838 | 33 Fr/35 Fr |
| | | 391-839 | 37 Fr/39 Fr |
| 391-833 | 13 Fr/15 Fr | 391-840 | 41 Fr/43 Fr |
| 391-834 | 17 Fr/19 Fr | 391-841 | 45 Fr/47 Fr |
| 391-835 | 21 Fr/23 Fr | 391-842 | 49 Fr/51 Fr |
| 391-836 | 25 Fr/27 Fr | 391-843 | 53 Fr/55 Fr |
| 391-837 | 29 Fr/31 Fr | 391-844 | 57 Fr/59 Fr |



Hegar Dilator Set €€

Stainless steel Hegar dilators provide the option of double-ended dilators. Available individually or in a set of 8.

| NUMBER | SIZE MM |
|--------|--|
| 61207 | 9 Fr (3 mm)/12 Fr (4 mm); 51 Fr (17 mm)/ 54 Fr (18 mm) (8/Set) |
| 61208 | 9 Fr (3 mm)/12 Fr (4 mm) |
| 61209 | 15 Fr (5 mm)/18 Fr (6 mm) |
| 61210 | 21 Fr (7 mm)/24 Fr (8 mm) |
| 61211 | 27 Fr (9 mm)/30 Fr (10 mm) |
| 61212 | 33 Fr (11 mm)/36 Fr (12 mm) |
| 61213 | 39 Fr (13 mm)/42 Fr (14 mm) |
| 61214 | 45 Fr (15 mm)/28 Fr (16 mm) |
| 61215 | 51 Fr (17 mm)/54 Fr (18 mm) |



Hank Cervical Dilators €€


Hank cervical dilators available individually or in a set of six ranging from 9 French (3 mm) to 20 French (6.66 mm). The set comes in a custom-fitted and numbered pouch.

| NUMBER | SIZE MM |
|--------|---------------------------------|
| 64-370 | 6 Dilators and Storage Pouch |
| 61201 | 9 Fr (3 mm)/10 Fr (3.33 mm) |
| 61202 | 11 Fr (3.67 mm)/12 Fr (4 mm) |
| 61203 | 13 Fr (4.33 mm)/14 Fr (4.67 mm) |
| 61204 | 15 Fr (5 mm)/16 Fr (5.33 mm) |
| 61205 | 17 Fr (5.67 mm)/18 Fr (6 mm) |
| 61206 | 19 Fr (6.33 mm)/20 Fr (6.67 mm) |



Mini Cervical Dilator Set €€

| NUMBER | DESCRIPTION |
|--------|--------------------------------------|
| 64-374 | 5 Dilators, Handle and Storage Pouch |



Small Cervical Dilator Set €€

| NUMBER | DESCRIPTION |
|--------|------------------------------|
| 64-372 | 6 Dilators and Storage Pouch |

Euro-Med™ Stainless Steel Endocervical Curettes



Kevorkian with Basket

| NUMBER | LENGTH |
|--------|------------|
| 66-442 | 12"/305 mm |



Handle Style

Kevorkian without Basket

| NUMBER | LENGTH |
|--------|------------|
| 66-440 | 12"/305 mm |



Kevorkian Collector™

| NUMBER | LENGTH |
|--------|------------|
| 66-441 | 12"/305 mm |



London™

| NUMBER | LENGTH |
|--------|--------------|
| 64-490 | 11.7"/297 mm |



Townsend™

| NUMBER | LENGTH |
|--------|--------------|
| 66-445 | 12.3"/312 mm |



Handle Style

Endocrette™

| NUMBER | LENGTH |
|--------|--------------|
| 66-447 | 12.3"/312 mm |

Euro-Med™ Endometrial Curettes



Handle Style

A.

B.



A. Novak

| NUMBER | SIZE | LENGTH |
|--------|------|-------------|
| 64-624 | 1 mm | 9.8"/248 mm |
| 64-626 | 2 mm | 9.8"/248 mm |
| 64-628 | 3 mm | 9.8"/248 mm |
| 64-630 | 4 mm | 9.8"/248 mm |

B. Randall

| NUMBER | SIZE | LENGTH |
|--------|--------|-------------|
| 64-620 | 3.9 mm | 9.9"/252 mm |



C.

Handle Style

C. Sims

| NUMBER | SIZE | LENGTH |
|---------|---------|------------|
| 392-110 | 6 mm | 11"/280 mm |
| 392-111 | 7 mm | 11"/280 mm |
| 392-112 | 8 mm | 11"/280 mm |
| 392-113 | 9 mm | 11"/280 mm |
| 392-114 | 11.5 mm | 11"/280 mm |
| 392-116 | 17 mm | 11"/280 mm |

Euro-Med™ Stainless Steel Hooks



Angle Hooks

| NUMBER | DESCRIPTION | LENGTH |
|--------|---------------|-------------|
| 36-630 | Straight Hook | 9.8"/249 mm |
| 36-632 | Right Hook | 9.6"/244 mm |



Iris Hook

| NUMBER | LENGTH |
|--------|--------------|
| 36-625 | 10.3"/262 mm |



Marlow Corkscrew Hook €

| NUMBER | LENGTH |
|--------|-------------|
| 67-451 | 9.3"/236 mm |



Double IUD Hook €

| NUMBER | LENGTH |
|---------|--------------|
| 392-151 | 10.6"/269 mm |



Emmett Tenaculum Hook €

| NUMBER | LENGTH |
|--------|-------------|
| 66-450 | 9.8"/249 mm |



IUD Hook €

| NUMBER | LENGTH |
|---------|--------------|
| 392-148 | 10.3"/262 mm |

Euro-Med™ Stainless Steel Cervical Tenacula



Braun Tenaculum Forceps €

| NUMBER | LENGTH |
|--------|-------------|
| 64-420 | 9.9"/252 mm |



Schroeder Tenaculum €

| NUMBER | LENGTH |
|--------|-------------|
| 64-410 | 9.9"/252 mm |



Duplay Tenaculum Double Curved €

| NUMBER | LENGTH |
|--------|------------|
| 61352 | 11"/280 mm |



Goldstein Grasp Cervical Stabilizer €

| NUMBER | LENGTH |
|---------|-------------|
| 64-6026 | 9.4"/239 mm |

Euro-Med™ Stainless Steel Dressing and Tissue Forceps



Allis €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61981 | 4 x 5 Teeth | 6" |
| 61983 | 5 x 6 Teeth | 7.5" |
| 61882 | 5 x 6 Teeth | 10" |



Potts-Smith Tissue €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61871 | 1 x 2 Teeth | 7" |
| 61872 | 1 x 2 Teeth | 8.2" |



Dressing Serrated €€

| NUMBER | LENGTH |
|--------|--------|
| 61804 | 5.6" |
| 61805 | 6.2" |



Potts-Smith Dressing €€

| NUMBER | LENGTH |
|--------|--------|
| 10-202 | 8.25" |



Adson Dressing Serrated €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|---------------|--------|
| 61830 | Serrated Tips | 4.75" |



Adson Tissue 1 x 2 Serrated €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|---------------------|--------|
| 61864 | 1 x 2 Serrated Tips | 4.75" |



Tissue 1 x 2 €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61841 | 1 x 2 Teeth | 5.75" |
| 61843 | 1 x 2 Teeth | 7" |
| 61844 | 1 x 2 Teeth | 8" |



Adson Tissue 1 x 2 €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61863 | 1 x 2 Teeth | 4.75" |



Adson-Brown Grasping €€

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 61866 | 7 x 7 Teeth | 4.75" |

Euro-Med™ Stainless Steel Uterine Forceps



Bozeman Uterine Forceps

| NUMBER | LENGTH |
|---------|--------------|
| 392-660 | 10.2"/257 mm |



Kelly Uterine Forceps

| NUMBER | LENGTH |
|---------|--------------|
| 417-131 | 12.5"/318 mm |



Foerster Uterine Polyp Forceps

| NUMBER | DESCRIPTION | LENGTH |
|---------|-------------|------------|
| 415-140 | Straight | 10"/254 mm |
| 415-141 | Curved | 10"/254 mm |

Euro-Med™ Stainless Steel Scissors



Operating Scissors (Straight)

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 62002 | B/B | 5.5" |
| 62008 | S/B | 5.5" |
| 62014 | S/S | 5.5" |
| 62004 | B/B | 6.5" |
| 62005 | B/B | 7.5" |
| 62010 | S/B | 6.5" |
| 62011 | S/B | 7.23" |



Spencer Suture Scissors

| NUMBER | DESCRIPTION | LENGTH |
|--------|-------------|--------|
| 62164 | Straight | 3.5" |

Euro-Med™ Stainless Steel Needle Holders



Mayo-Hegar Needle Holder (Wider Jaw)

| NUMBER | LENGTH |
|--------|--------|
| 62809 | 5.5" |
| 62810 | 6.3" |
| 62815 | 7.1" |
| 62820 | 8.1" |



Heaney Needle Holder (Curved)

| NUMBER | LENGTH |
|--------|--------|
| 62825 | 8.3" |



Crile-Wood Needle Holder (Narrow Jaw)

| NUMBER | LENGTH |
|--------|--------|
| 62840 | 7.1" |



Metzenbaum Scissors (Straight or Curved)

| Straight | | Curved | |
|----------|--------|--------|--------|
| NUMBER | LENGTH | NUMBER | LENGTH |
| 62116 | 7" | 62131 | 5.9" |
| 62127 | 8" | 62118 | 7" |
| 62120 | 9.1" | 62119 | 8.1" |
| 62122 | 11.2" | 62124 | 9" |
| | | 62126 | 11" |



Mayo Scissors (Straight or Curved)

| Straight | | Curved | |
|----------|--------|--------|--------|
| NUMBER | LENGTH | NUMBER | LENGTH |
| 62063 | 5.7" | 62067 | 5.5" |
| 62065 | 6.8" | 62069 | 6.75" |
| 62066 | 9" | 62070 | 9" |



Iris Scissors (Straight or Curved)

| NUMBER | DESCRIPTION | LENGTH |
|--------|--------------|--------|
| 62140 | Straight S/S | 4.5" |
| 62143 | Curved S/S | 4.5" |

Endometrial Samplers



Pipelle® Endometrial Suction Curette CA* CE

| NUMBER | DESCRIPTION |
|--------|--|
| 8200 | Pipelle Endometrial Suction Curette (25/Box) |



Pipet Curet™ Endometrial Suction Curette CA* CE

| NUMBER | DESCRIPTION |
|--------|--|
| MX140 | 3 mm (50/Box) – Greater integrity for stenotic Os |
| MX145 | 3 mm Soft Yellow (50/Box) – Pliable for molding catheter |
| MX150 | 4 mm (50/Box) – Multiparous patient |



Uterine Explora™ Model I Curette with Sharp Randall-Type Cutting Edge and Vacu-Lok Syringe CA* CE

Disposable Explora Model incorporates a long nylon curette with Randall-type cutting edge, stylet and Vacu-Lok syringe that adjusts suction from approximately 15- to 20-inch Hg negative pressure (6 cc to 12 cc).

| NUMBER | DESCRIPTION |
|--------|---|
| MX120 | Uterine Explora with Vacu-Lok Syringe, 3 mm, Model I (25/Box) |
| MX124 | Uterine Explora with Vacu-Lok Syringe, 4 mm, Model I (25/Box) |



Uterine Explora™ Model II Curette with Dual Hole Design and Vacu-Lok Syringe CA* CE

Uterine Explora Model II offers the same features as the Model I plus it has a dual-port opening. The Vacu-Lok syringe allows you to adjust the suction from approximately 15- to 20-inch Hg negative pressure (6 cc to 12 cc).

| NUMBER | DESCRIPTION |
|--------|--|
| MX130 | Uterine Explora II with Vacu-Lok Syringe Model II (25/Box) |




Disposable Uterine Sound and Dilators



Fundus Sound™ XM CA☼

| NUMBER | DESCRIPTION |
|--------|---|
| 1196 | Fundus Sound XM Disposable Uterine Sound (25/Box) |



Milex® Dilator / Sound Set CA☼

| NUMBER | DESCRIPTION |
|--------|--|
| MX11 | Milex Os Locator (25/Box) |
| MX12 | Milex Canal Dilator (25/Box) |
| MX13 | Milex Fundal Sound (25/Box) |
| MX21 | Milex Dilator/Sound Set, 3 dilators per set: MX11, MX12, MX13 (10/Box) |



Os Finder™ CA☼ CE


| NUMBER | DESCRIPTION |
|--------|--------------------|
| 1176 | Os Finder (50/Box) |

Dilators

Comfort-Flex Dilator® Set CA*

- Os Locator (1 ring) – opens the external Os
- Canal Dilator (2 rings) – traces the natural path of the cervical canal while steadily dilating
- Fundus Sound™ (3 rings) – opens the internal Os and sounds the fundus to determine uterine depth

| NUMBER | DESCRIPTION |
|--------|--|
| 64-376 | Comfort-Flex Dilator Set, 3 dilators per set (Os Locator, Canal Dilator, Fundus Sound) |



Milex® Vaginal-Hymenal Silicone Dilators CA*

| NUMBER | DESCRIPTION |
|---------|--|
| MX20 | Milex Vaginal Dilator – Set of 4: Small, Medium, Large, Extra-Large with Trimo-San (USA sales only) |
| MX20-CA | Milex Vaginal Dilator – Set of 4: Small, Medium, Large, Extra-Large without Trimo-San |
| MX20I | Milex Vaginal Dilator Single – Small 5 ¹ / ₈ " x 1 ⁵ / ₁₆ " |
| MX20IM | Milex Vaginal Dilator Single – Medium 5 ¹ / ₂ " x 1 ¹ / ₁₆ " |
| MX20IL | Milex Vaginal Dilator Single – Large 5 ³ / ₄ " x 1 ¹ / ₄ " |
| MX20IXL | Milex Vaginal Dilator Single – Extra-Large 6 ¹ / ₈ " x 1 ¹ / ₂ " |

Milex® Vaginal-Rectal Silicone Dilators CA*

| NUMBER | DESCRIPTION |
|---------|---|
| MX25 | Milex Rectal Dilator – Set of 4: Small Rectal, Medium Rectal, Large Rectal, Extra-Large Vaginal-Rectal with Trimo-San (USA sales only) |
| MX25-CA | Milex Rectal Dilator – Set of 4: Small Rectal, Medium Rectal, Large Rectal, Extra-Large Vaginal-Rectal without Trimo-San (International sales only) |
| MX25IS | Milex Vaginal-Rectal Dilator Single – Small 3" x 1 ¹ / ₁₆ " |
| MX25IM | Milex Rectal Dilator Single – Medium 3 ¹ / ₈ " x 1 ⁵ / ₁₆ " |
| MX25IL | Milex Rectal Dilator Single – Large 3 ⁹ / ₁₆ " x 1 ¹ / ₈ " |
| MX25IXL | Milex Vaginal-Rectal Dilator Single – Extra-Large 3 ¹⁵ / ₁₆ " x 1 ⁵ / ₁₆ " |

Milex® Vaginal Silicone Dilators CA*

| NUMBER | DESCRIPTION |
|---------|--|
| MX23 | Milex Dilator – Set of 4: Small, Medium, Large, Extra-Large with Trimo-San (USA sales only) |
| MX23-CA | Milex Dilator – Set of 4: Small, Medium, Large, Extra-Large without Trimo-San |
| MX23I | Milex Dilator Single – Small 2 ¹ / ₁₆ " x 1 ¹ / ₂ " x 3 ³ / ₈ " |
| MX23IM | Milex Dilator Single – Medium 2 ⁵ / ₁₆ " x 9 ⁹ / ₁₆ " x 7 ⁷ / ₁₆ " |
| MX23IL | Milex Dilator Single – Large 2 ¹¹ / ₁₆ " x 5 ⁵ / ₈ " x 1 ¹ / ₂ " |
| MX23IXL | Milex Dilator Single – Extra-Large 3" x 1 ¹¹ / ₁₆ " x 9 ⁹ / ₁₆ " |

Cervical and Endocervical Samplers

Medscand® Cytobrush Plus® Kits

| NUMBER | DESCRIPTION |
|-------------|--|
| 10115-10/50 | Cytobrush Plus Kits – Cell Collector (sets include brush, scraper, fixative, slide and slide transporter); 10 Trays with 50 Pap Smear test sets (500 total Tests/Case) |
| 60115-10/50 | Cytobrush Plus Cell Collector – Gentle Touch, 10 Trays with 50 Pap Smear test sets (500 total Tests/Case) |
| 60115-20/50 | Cytobrush Plus Cell Collector – Gentle Touch, 20 Trays with 50 Pap Smear test sets (1,000 total Tests/Case) |



Medscand® Cytobrush Plus® CA* CE

| NUMBER | DESCRIPTION |
|--------|--|
| C0021 | Cytobrush Plus – 2 Bags, 500 Brushes each (1,000/Box) |
| C0004 | Cytobrush Plus – 10 Bags, 10 Brushes each (100/Box) |
| C0005 | Cytobrush Plus – 25 Bags, 100 Brushes each (2,500/Box) |
| C0012 | Cytobrush Plus Sterile, 1 Brush per Pouch (40/Box) |

Medscand® Cytobrush Plus® GT CA* CE

| NUMBER | DESCRIPTION |
|--------|--|
| C0121 | Cytobrush Plus GT – 2 Bags, 500 Brushes each (1,000/Box) |
| C0104 | Cytobrush Plus GT – 10 Bags, 10 Brushes each (100/Box) |
| C0105 | Cytobrush Plus GT – 25 Bags, 100 Brushes each (2,500/Box) |
| C0112 | Cytobrush Plus GT Sterile, 1 Brush per Pouch (40/Box) |
| C0305 | Cytobrush Plus GT Scored – 25 Bags, 100 Brushes each (2,500/Box) |



Medscand® Pap Perfect® Spatula

| NUMBER | DESCRIPTION |
|--------|--|
| 1108 | Pap-Perfect Endocervical Scrapers (Wooden) (500/Box) |
| 11080 | Pap-Perfect Plastic Spatulas (500/Box) CA* CE |

Diagnostic Disposables



Papette® Cervical Cell Collector

| NUMBER | DESCRIPTION | |
|--------|--|--------|
| 908001 | Papette Cervical Cell Sampler Brush (100 per Bag) | CA☼ CE |
| 908003 | Papette Cervical Cell Sampler Brush (10 per Box; 1,000 pieces) | CA☼ CE |
| 908004 | Papette Cervical Cell Sampler Brush (40 per Box; 4,000 pieces) | CA☼ CE |

Endocell® Disposable Endometrial Cell Sampler

| NUMBER | DESCRIPTION | |
|----------|--|--------|
| 908014 A | Endocell – 35 per box | CA☼ CE |
| 908015 | Endocell – 6 boxes of 35 (210 pieces) | CA☼ CE |
| 908016 | Endocell – 15 boxes of 35 (525 pieces) | CA☼ CE |



The Inseminator

Intrauterine Insemination Device for Assisted Reproduction

- 3.5 French polypropylene catheter enlarged in the proximal area to accommodate the syringe tip
- Comes in 4.5 in. or 5.5 in. length
- 4.5 in. may be purchased without syringe

| NUMBER | DESCRIPTION | |
|----------|--|-----|
| 920012 | The Inseminator 5.5 in. (50 per Box) | |
| 920012-1 | The Inseminator 4.5 in. (50 per Box) | CA☼ |
| 920012-2 | The Inseminator 4.5 in. without syringe (50 per Box) | CA☼ |



Flexisound™ Intrauterine Sound

Graduated in centimeters.

| NUMBER | DESCRIPTION |
|--------|--|
| 920008 | Flexisound Disposable Intrauterine Sounds (25 per Box) |



Disposable Os Finders

Disposable Os Finders

- Os Finder
- Canal Finder
- Fundus Finder

| NUMBER | DESCRIPTION |
|--------|--|
| 920019 | Disposable Os Finders – Set of 3 (10 Sets per Box) |



Breast Aspiration Dual-Port Biopsy Needle

| NUMBER | DESCRIPTION | |
|--------|---|-----|
| MX40 | Dual-Port Needle with 3-finger control with 12 cc Syringe, 3" (12/Box) | |
| MX41 | Dual-Port Needle with 3-finger control with 12 cc Syringe, 1½" (12/Box) | CA♣ |
| MX42 | Dual-Port Needle, No Syringe, 1½" (5/Box) | CA♣ |



Word Catheter CA♣


| NUMBER | DESCRIPTION |
|--------|--|
| 53404 | Word Bartholin Gland Catheter. Includes syringe and scalpel. Latex-free. (6/Box) |



Milex® Cervical Spatulas

| NUMBER | DESCRIPTION |
|--------|--|
| MX600 | Milex Cervical Spatulas, Small, Yellow (250/Box) |
| MX610 | Milex Cervical Spatulas, Medium, White (250/Box) |

Leisegang™ LM-900 Cryosurgery System



| NUMBER | DESCRIPTION | |
|--------|--|---|
| 50501 | LM-900 N ₂ O Cryosurgery Unit with Case | € |
| 50502 | LM-900 CO ₂ Cryosurgery Unit with Case | € |

Accessories

| NUMBER | DESCRIPTION | |
|--------|-------------------------------------|---|
| 50000 | Empty 20 lb. N ₂ O | |
| 50010 | Empty 20 lb. CO ₂ | |
| 50003 | 20 lb. Cylinder Cart | |
| 50007 | E-Cylinder Adapter N ₂ O | € |
| 50008 | E-Cylinder Adapter CO ₂ | € |
| 50009 | Carrying Case | |
| 50011 | O-Rings (10 μ) | € |
| 50012 | Cleaning Plugs (1 μ) | € |

LEISEGANG LM-900 Probes CA* CE

| NUMBER | DESCRIPTION |
|--------|--|
| 50101 | Microderm - 2 mm D x 25 mm L |
| 50102 | Soll Lash Epilating - 2 mm D x 27 mm L |
| 50103 | Condyloma - 5 mm D x 25 mm L |
| 50104 | Derm Bevel - 5 mm D x 24 mm L |
| 50105 | General Purpose - 19 mm D x 8 mm L |
| 50106 | Large Derm Bevel |
| 50107 | Curved Anorectal - 8 mm D x 38 mm L |
| 50108 | General Purpose - 19 mm D x 8 mm L |

| NUMBER | DESCRIPTION |
|--------|---------------------------------------|
| 50109 | Exocervical - 19 mm D x 5 mm L |
| 50110 | Exo/Endocervical - 19 mm D x 10 mm L |
| 50111 | Special Exo/Endo - 19 mm D x 10 mm L |
| 50112 | Exo/Endocervical - 19 mm D x 20 mm L |
| 50113 | Maxiderm - 25 mm D x 8 mm L |
| 50114 | Large Exocervical - 25 mm D x 10 mm L |
| 50115 | Jumbo Exo/Endo - 25 mm D x 17 mm L |

Cryosurgical Systems:

One- and Two-Trigger Systems Provide Conductive, Reusable Cryotips and Disposable Plastic Shields Available for Cryosurgical Treatment



LL100™ System (N₂O) Two-trigger

Freeze and defrost only when triggered. Use with cryotips and disposable plastic shields.

LL100™ System (N₂O) & LLCO₂™ System

| NUMBER | DESCRIPTION | |
|--------|--|--------|
| 900001 | LL100 Cryosurgical System – N ₂ O (Freezer Only) | |
| 900160 | LLCO ₂ Cryosurgical System – CO ₂ (Freezer Only) | |
| 901060 | N ₂ O Cylinder (unfilled) | CA* CE |
| 901061 | CO ₂ Cylinder (unfilled) | CE |
| 901011 | Cart for 20 lb. Cylinder | |

LLCO₂™ System Single-trigger for CO₂ Gas

- Squeeze trigger to freeze; release to defrost
- Use with medical-grade CO₂ gas for optimal performance
- Use with CooperSurgical's standard cryotips with disposable plastic shields

Both the LL100 and the LLCO₂ are available in three-tip package specials for additional savings. Choose three cryotips from our comprehensive offerings when ordering. The following are recommended for general/family practice procedures: T-1905, T-2514, and T-0500 for GYN procedures and the T-0524, T-0219, and T-1900.



WA1000B™ Cryo Console System

Uses LL100™ two-trigger or LLCO₂™ single-trigger freezers

- Temperature Indicator with integrated timer
- Cart carries a single 20 lb. cylinder
- Use with CooperSurgical's large selection of cryotips and shields

| NUMBER | DESCRIPTION | |
|--------------|---|--------|
| 900506-1-N2O | WA1000B Cryosurgical Console – Blank N ₂ O LL100 | CA* CE |
| 900506-1-CO2 | WA1000B Cryosurgical Console – Blank CO ₂ LL100 | CA* |
| 900506-2-N2O | WA1000B Cryosurgical Console – Digital Temp. (N ₂ O) Indicator, Timer, Base, and LL100 | CE |
| 900506-2-CO2 | WA1000B Cryosurgical Console – Digital Temp. (CO ₂) Indicator, Timer, Base, and LLCO ₂ | |




Cryosurgical Tips for Use with LL100™ and LLCO₂™ CA* CE

Conductive, thermal, reusable, sterilizable cryotips with disposable plastic shields.

| NUMBER | DESCRIPTION | NUMBER | DESCRIPTION | NUMBER | DESCRIPTION | NUMBER | DESCRIPTION |
|----------|-------------|----------|-------------|----------|-------------|----------|---|
| 900300AA | (T-0200) | 900200AA | (T-0219) | 900206AA | (T-0826) | 900212AA | (T-1920) |
| 900301AA | (T-0500) | 900201AA | (T-0219) | 900207AA | (T-0823) | 900213AA | (T-2500) |
| 900302AA | (T-0700) | 900202AA | (T-0524) | 900208AA | (T-1900) | 900214AA | (T-2507) |
| 900303AA | (T-0900) | 900203AA | (T-0519) | 900209AA | (T-1905) | 900215AA | (T-2514) |
| 900305AA | (T-1300) | 900204AA | (T-0832) | 900210AA | (T-1920) | 900217AA | (T-0832) |
| 900306AA | (T-1500) | 900205AA | (T-0819) | 900211AA | (T-1910) | 900221AA | Disposable Plastic Shields for Cryotips (100 per Bag) |

ABBI®

Air Bubble Based Infuser



ABBI® Air Bubble Based Infuser

| NUMBER | DESCRIPTION |
|----------|---|
| ABBI | ABBI (Air Bubble Based Infuser) – Single-Use Disposable (1/Box) |
| ABBI-Kit | ABBI (Air Bubble Based Infuser) with H/S Elliptosphere Catheter (1/Box) |



H/S Elliptosphere™ Catheter Set and Tray

Sonohysterography and Hysterosalpingography

Available in two versions – with and without an integrated stylet – the H/S Elliptosphere features latex-free construction.

| NUMBER | DESCRIPTION | |
|----------|--|--------|
| 61-4005 | H/S Elliptosphere Catheter, 5 Fr (10/Box) | CA* CE |
| 61-4005S | H/S Elliptosphere Catheter Set with Integrated Stylet (10/Box) | CA* CE |
| 61-4205 | H/S Elliptosphere Tray with 5 Fr Catheter (10/Case) | |



61-4205
H/S Elliptosphere Tray

H/S Catheters


For Visualization of Uterine Pathology

H/S Catheter Set CA* CE

Designed for use in sonohysterography and hysterosalpingography, the catheter can be used with oil- or saline-based media depending on the procedure requirements.

- Available with integrated stylet
- Comes in both 5 and 7 French


| NUMBER | DESCRIPTION |
|---------|--|
| 61-5005 | H/S Catheter, Kraton, 5 Fr (10/Box) |
| 61-5007 | H/S Catheter, Kraton, 7 Fr (10/Box) |
| 61-0020 | H/S Stylet Catheter Stiffener (10/Box) |



H/S Procedure Tray

- Fenestrated drape
- Disposable speculum
- Disposable cervical dilator
- Three swab sticks
- Sterile wrap
- Inner tray for povidone-iodine solution
- One pack of povidone-iodine solution
- 36-in. extension tube
- 20 cc syringe with 18-gauge needle
- Lubricating jelly packet
- Sanitary napkin

| NUMBER | DESCRIPTION |
|---------|---|
| 61-5205 | H/S Procedure Tray – Round 5 Fr |
| 61-5207 | H/S Procedure Tray – Round 7 Fr |
| 61-3205 | H/S Procedure Tray – Round 5 Fr with Stylet |
| 61-3207 | H/S Procedure Tray – Round 7 Fr with Stylet |



Tampa Infusion Catheter CA* CE

For sonohysterography – ultrasound examination of the uterus.

| NUMBER | DESCRIPTION |
|---------|-------------------------|
| 61-2005 | Tampa Catheter (20/Box) |



HSG Procedure Tray

The HSG Procedure Tray is for use in hysterosalpingography and sonohysterography procedures.

| NUMBER | DESCRIPTION |
|--------|----------------|
| 6050T | Procedure Tray |

Fertility Assessment



SelectMucus® Endocervical Aspirator

| NUMBER | DESCRIPTION | |
|--------|---|-----|
| SM301 | SelectMucus Endocervical Sampler (25/Box) | CA* |



Aspirette® Endocervical Aspirator CA* Cervical Sampling Made Easy

With a 3.1 mm O.D. body and rounded distal tip, cervical insertion is easy.

| NUMBER | DESCRIPTION |
|--------|---|
| 8700 | Aspirette Endocervical Aspirator (25/Box) |



 **Paragard**[®]
intrauterine copper contraceptive

by CooperSurgical[®]

Ordering Paragard

Two simple ways to order the Paragard IUS for your patients:



Paragard DirectSM is our “buy and bill” option for ordering Paragard that allows you to take advantage of volume discount pricing.



Create an account



Order online, by fax or by phone



Ship to one or more locations

Visit ParagardDirect.com or call **1-877-PARAGARD** to learn more.



With the Paragard Specialty PharmacySM program, your patient’s prescription for Paragard is fulfilled through one of our specialty pharmacies and shipped directly to your office.



Process the claim with your patient’s health insurer



Collect any applicable co-pay for the unit, if needed



Ship Paragard directly to your office

To order through a specialty pharmacy:

1. Visit ParagardAccessCenter.com to download our specialty pharmacy request form
2. Complete the form and fax it back to the selected specialty pharmacy



Designed for the collection of cord blood and cord tissue stem cells.

Ordering a CBR® Newborn Stem Cell Preservation Collection Kit

Have your patient visit cordblood.com or call CBR® (Cord Blood RegistrySM) directly at 1.888.267.3256 or you can email providersupport@cordblood.com to request collection kits for your office or hospital.

Fetal Monitor



F9 Fetal Monitor

Fetal monitor with 2 ultrasound transducers, 1 toco transducer, 1 event marker, 2 reusable belts, ultrasound gel, printer paper, power cord, system user manual.
Optional: CooperSurgical Insight™ software, RS232 and ethernet cables, rechargeable lithium-ion battery.

| NUMBER | DESCRIPTION |
|-----------------|--|
| MS9-108022-UL-R | CooperSurgical Fetal Monitor F9 |
| MS9-108022-SWP | CooperSurgical Insight Software Package (Includes software and cables) |
| M21-64072 | Rechargeable Lithium-ion Battery |
| MS9-107634 | Deluxe Cart |
| MS9-108166 | Base Cart with Basket and Laptop Tray |
| MS9-108108 | Base Cart with Basket |
| MS3-31528-A1 | Transducer 1 (pink, 4-pin) |
| MS3-31528-A1-R | Transducer 1 (pink, 6-pin) |
| MS3-107705-A1 | Transducer 2 (yellow, 4-pin) |
| MS3-107705-A1-R | Transducer 2 (yellow, 6-pin) |
| MS3-31527-A1 | TOCO Probe |
| MS3-31112 | Event Marker |
| MS1R-02264 | Transducer Belt (reusable) |
| M25R-75111 | Recording Paper (90 mm x152 mm x 150 p, US standard) (Pack of 3) |



Fetal2EMR™

Fetal Monitor with the capability of transferring NSTs to electronic medical records.

| NUMBER | DESCRIPTION |
|----------|---|
| 902300-R | Fetal2EMR Monitor and Accessories |
| 902320 | Fetal2EMR Insight™ Software; Includes Installation CD and Connecting Cables |
| 902302 | FHR/Ultrasound Transducer (4-pin) |
| 902302-R | FHR/Ultrasound Transducer (6-pin) |
| 902304 | Thermo Paper 90 mm x 152 mm x 150 p (Pack of 3) |
| 902340 | Fetal2EMR Rolling Cart with Tray |
| 902341 | Fetal2EMR Rolling Cart without Tray |

Obstetrical Dopplers and Probes



LifeDop® 150

Handheld non-display Doppler with rechargeable battery and audio recording options. Compatible with all LifeDop non-directional probes.

| NUMBER | DESCRIPTION | |
|----------|---|--------|
| L150-SDW | LifeDop 150 with Waterproof 2 MHz Probe | CA* CE |

LifeDop® 250 Series

Handheld Dopplers with large displays for indication of heart rate and battery status. Rechargeable batteries and audio recording options are available on specific models. Compatible with all LifeDop non-directional probes.

| NUMBER | DESCRIPTION | |
|------------|---|--------|
| L250-SD2 | LifeDop 250 Disp. with 2 MHz Probe | CA* CE |
| L250-SD3 | LifeDop 250 Disp. with 3 MHz Probe | CA* CE |
| L250-SDW | LifeDop 250 Disp. with 2 MHz Waterproof Probe | CA* CE |
| L250AR-SDW | LifeDop 250 Disp. with Audio Rec., Recharge & 2MHz Waterproof | CA* CE |
| L250R-SD2 | LifeDop 250 Disp. with Recharge & 2 MHz Probe | CA* CE |
| L250R-SD3 | LifeDop 250 Disp. with Recharge & 3 MHz Probe | CA* CE |
| L250R-SDW | LifeDop 250 with Recharge & 2 MHz Waterproof Probe | |
| K200 | Roll Stand with Storage Basket for L150 or L250 | |



LifeDop® 350 Series

The LifeDop 350 Doppler is a tabletop Doppler system designed to detect the fetal heart rate and monitor blood flow. It has a large heart rate display, comes with rechargeable batteries, and accepts all LifeDop non-directional probes.

- Two probe holders for multiple applications
- User-accessible batteries
- Stand and wall mounting available

| NUMBER | DESCRIPTION | |
|-----------|--|--------|
| L350R-SD2 | LifeDop Tabletop Doppler with 2 MHz Probe | CA* CE |
| L350R-SD3 | LifeDop Tabletop Doppler with 3 MHz Probe | CA* CE |
| L350R-SDW | LifeDop Tabletop Doppler with 2 MHz Waterproof Probe | CA* CE |
| K310 | Wall Mount for L350R | |
| K220 | Roll Stand with Storage Basket for L350R | |



LifeDop® Non-Directional Obstetrical Probes

All LifeDop non-directional probes are completely interchangeable for use on the L150 and the L250 and L350 series of Dopplers. Most models have user-replaceable cables.

2 MHz OB Probe Wide-angle beam for deep penetration needed in late term pregnancy and larger patients.

3 MHz OB Probe The general-purpose obstetrical probe for fetal heart detection, from earliest stages through delivery.

2 MHz Waterproof OB Probe Waterproof design for water-assisted labor and delivery applications.

| NUMBER | DESCRIPTION | |
|--------|--|--------|
| SD2 | 2 MHz OB Probe | CA* CE |
| SD3 | 3 MHz OB Probe | CA* CE |
| SDW | 2 MHz Waterproof OB Probe | CA* CE |
| C250CC | Obstetrical Vascular Combination – L250 Doppler, 3 MHz and 8 MHz probe | CA* CE |
| C250CB | Obstetrical Combination – L250 Doppler, 2 MHz and 3 MHz probe | CA* CE |

Vascular Assessment Dopplers and Probes



LifeDop® 150 Series

These handheld, non-display Dopplers for vascular assessment offer full line of non-directional LifeDop probes. Rechargeable batteries and audio recording options are available on certain models.

| NUMBER | DESCRIPTION |
|-----------|---|
| L150R-SD4 | LifeDop 150 with 4 MHz Probe, Rechargeable |
| L150R-SD5 | LifeDop 150 with 5 MHz Probe, Rechargeable |
| L150R-SD8 | LifeDop 150 with 8 MHz Probe, Rechargeable |
| L150R-STN | LifeDop 150 with 8 MHz Sterilizable Probe, Rechargeable |
| STNPK1 | L150 Doppler with Sterilizable Probe and Stand |

LifeDop® 250

The handheld LifeDop L250 Doppler has a large display for heart rate and battery status.

| NUMBER | DESCRIPTION |
|----------|---|
| L250-STN | LifeDop 250 Disp. with 8 MHz Sterilizable Probe |



LifeDop® 350

The LifeDop 350 has a large heart rate display, rechargeable batteries, and accepts all LifeDop non-directional probes. Enhanced features include:

- Two probe holders for multiple applications
- User accessible batteries
- Wall mounting

| NUMBER | DESCRIPTION |
|--------|---|
| K310 | Wall Mount for Added Convenience and Security |



LifeDop® Vascular Assessment Probes

CA★ CE

LifeDop probes are interchangeable for use on the L150 Series and L350 Dopplers, with the exception of the sterilizable probes, which have some restrictions. Most models have user replaceable cables.

- 4 MHz Vascular Probe
- 5 MHz Vascular Probe
- 8 MHz Vascular Probe
- 8 MHz Sterilizable Probe

| NUMBER | DESCRIPTION |
|--------|-----------------------------------|
| SD4 | 4 MHz Vascular Probe |
| SD5 | 5 MHz Vascular Probe |
| SD8 | 8 MHz Vascular Probe |
| STN | 8 MHz Sterilizable Vascular Probe |

Full-Featured ABI Systems for Peripheral Arterial Disease Diagnosis

Vista AVS®

A Full-Featured ABI/Vascular System

The Vista AVS offers three modalities for obtaining systolic pressures and waveforms and is designed to accommodate both standard and custom protocols. The accompanying compact printer prints two ankle waveforms along with the ABI results. The system includes:

- 3 modalities: PVR, bi-directional Doppler probe and PPG probe
- Automatic ABI, TBI and segmental calculation
- Performs seated ABI on mobility impaired patients
- Software for full-page reporting, exam storage and interface with EMR
- 9 blood pressure cuffs: (2) 10 cm, (4) 12 cm, (2) 17 cm contoured thigh cuffs, (1) digit cuff
- Automatic inflation/deflation
- Graphic waveform display with step-by-step ABI instructions

Optional Features for the Vista AVS®

5 MHz Bi-directional Probe

For deeper vessel penetration on larger patients.

Cuff Selector

Up to 10 vascular cuffs can be connected to this unit at the same time, for upper and lower extremity or digit pressure measurements.

| NUMBER | DESCRIPTION | |
|----------|---|-----|
| L500VA | Full-featured ABI/Vascular System | CA☐ |
| SD5V | 5 MHz Bi-directional Probe | |
| K290 | Cuff Selector | |
| G150 | Ultrasound Gel (60 grams) | |
| VISTCAP1 | Extended one year warranty for Vista system | |



Vista ABI™

This system includes:

- 2 modalities: PVR and bi-directional Doppler probe
- Automated cuff inflation/deflation
- Built-in printer
- Graphic waveform display with step-by-step ABI instructions
- 4 cuffs: (2) 10 cm, (2) 12 cm

Optional Features for the Vista ABI™

PPG Package

A sensitive arterial PPG probe for obtaining both toe and limb pressures. Comes with Velcro to secure placement on the digit. This package includes a PPG probe and a 1.9 cm digit cuff for performing TBI exams on patients with non-compressible arteries.

PC Download Software

An optional software package for the Vista ABI with full-page reporting and exam storage includes download capability for integration with EMR systems (.pdf and .csv file formats – USB Interface to PC).

| NUMBER | DESCRIPTION | |
|--------|----------------------------|-----|
| L450VA | Vista ABI System | CA☐ |
| L450PG | PPG Package | |
| L450SF | PC Download Software | |
| SD5V | 5 MHz Bi-directional Probe | |
| G150 | Ultrasound Gel (60 grams) | |

*The Vista ABI does not accommodate segmental exams.



Portable and Handheld ABI Systems



LifeDop® 300 ABI System with Four Cuffs and Aneroid

A Handheld ABI System with ABI Calculation and
Dual Waveform Printout

Portable bi-directional Doppler with keypad for entering pressure
data and calculating the indices. The accompanying compact
printer prints two ankle waveforms along with the ABI results.

System includes:

- Printer
- 8 MHz bi-directional probe
- Adhesive-backed label paper
- ABI chart, forms
- ABI video, examination guide, quick
reference guide, and ABI glossary
- Carrying case

| NUMBER | DESCRIPTION | |
|--------|---|-------|
| L300AC | LifeDop 300 ABI System with 4 Cuffs and Aneroid | CA* € |
| L250AC | LifeDop 250 ABI System with 4 Cuffs and Aneroid | CA* € |
| L250AB | LifeDop 250 ABI System with 2 Cuffs and Aneroid | CA* € |
| SD5B | 5 MHz Bi-directional Probe | |
| SD8B | 8 MHz Bi-directional Probe | |
| K180 | Adhesive-backed label paper | |
| G150 | Ultrasound Gel | |



LifeDop® 250 ABI System with Cuffs and Aneroid

Portable bi-directional Doppler with compact
printer. The system includes handheld Doppler for
diagnosing PAD, cuffs and aneroid package, printer
to capture waveforms, an 8 MHz bi-directional
probe, adhesive-backed label paper, carrying case,
and ABI reference documents and video.



VANMAX Shown

Vantage ABI™

The Vantage ABI™ is a cuff-based system that uses Digital Fourier
Oscillometric or DFO waveform analysis. The accompanying compact
printer prints two ankle waveforms along with the ABI results.

VANABI – Vantage ABI Package 1

- Vantage ABI
- Graphic display with convenient touch-screen to navigate the system
- USB memory stick for data transfer and storage (EMR interface)
- Four blood pressure cuffs: (2) 10 cm cuffs for ankles R/L and (2)
12 cm cuffs for arms R/L
- Automatic inflation/deflation

VANMAX – Vantage ABI Package 2

- Vantage ABI
- Graphic display with touch-screen to navigate the system
- USB memory stick for data transfer and storage (EMR interface)
- Built-in printer for printing ankle waveforms on adhesive-backed
label paper
- Four blood pressure cuffs: (2) 10 cm cuffs for ankles R/L and
(2) 12 cm cuffs for arms R/L
- Automatic inflation/deflation

| NUMBER | DESCRIPTION |
|--------|--|
| VANABI | Vantage ABI – Package 1 |
| VANMAX | Vantage ABI with Stand – Package 2 |
| K180 | Adhesive-backed label paper |
| C019 | Cuff and Hose Assembly – 1.9 cm Digit Cuff (VISTA Systems Only) |
| C100 | Cuff and Hose Assembly – 10 cm |
| C101 | Cuff and Hose Assembly – 10 cm – Ankle Left |
| C102 | Cuff and Hose Assembly – 10 cm – Ankle Right |
| C120 | Cuff and Hose Assembly – 12 cm |
| C121 | Cuff and Hose Assembly – 12 cm – Arm Left |
| C122 | Cuff and Hose Assembly – 12 cm – Arm Right |
| C170 | Cuff and Hose Assembly – 17 cm Contoured Thigh Cuff |

Ophthalmic Cryosurgery

Insulating Sleeves

| NO. | DESCRIPTION |
|-------|--|
| P805 | For Probes 124/125/142/161/162/166/167 (Retinal) (6/Pkg) |
| 805-1 | For Probes 143/144 (Hammerhead) (6/Pkg) |
| 32709 | For Probes 128/132/147 (Micro) (1/Pkg) |
| 32710 | For Probes 126/127 (Cataract) (6/Pkg) |
| 32787 | For Probes 139/163 (Retinal Spherical) (1/Pkg) |
| 33098 | For Probes 130/164/165 (Glaucoma) (6/Pkg) |
| 33555 | For Probes 120/121 (Retinal Pediatric) (1/Pkg) |



Wallach® UltraFreeze™

Portable cryosurgical system for most office skin procedures. The system includes an assortment of luer lock open spray apertures. A plastic reusable plate with six conical openings is available to help direct liquid nitrogen spray in concentrated areas. Available in .5 liter.

Wallach UltraFreeze™

| NUMBER | DESCRIPTION |
|----------|---|
| 900108 | Accu-Shield Cone Spray Plate (set of 6 varied cones) |
| 900111-1 | Manual Discharge Device for 10 liter Dewar |
| 900111-3 | Manual Discharge Device for 20 liter Dewar |
| 900112 | UltraFreeze Stand |
| 900076 | UltraFreeze Liquid Nitrogen Sprayer (.5 liter)* |
| 900113 | Aperture kit for UltraFreeze (4 spray apertures and one 20 gauge bent needle) |



Index

| | |
|---|-------|
| Biopsy Punches | |
| Euro-Med..... | 38-41 |
| Cervical Samplers | 65 |
| Colposcopy | |
| OptiK Model 1 | 31 |
| OptiK Model 2 | 32-33 |
| Wallach® Zoom Magnification Colposcopes..... | 34 |
| Cord Blood Preservation | |
| CBR | 80 |
| Cryosurgery | |
| Leisegang™ LM-900 Cryosurgery System..... | 70-71 |
| Cryosurgical Systems | 72-73 |
| Ophthalmic Cryosurgery | 90 |
| UltraFreze | 91 |
| Cytobrush | 65 |
| Curettes-Endocervical | |
| Stainless Steel | 52 |
| Curettes-Endometrial | |
| Disposable..... | 66 |
| Stainless Steel | 53 |
| Dilators | |
| Milex Vaginal-Hymenal Silicone Dilators..... | 64 |
| Milex Vaginal-Rectal Silicone Dilators..... | 64 |
| Milex Vaginal Silicone Dilators | 64 |
| Milex Dilator/Sound Set..... | 64 |
| Dopplers | |
| Obstetrical Doppler Probes and Accessories | 82-83 |
| Vascular Dopplers..... | 84-87 |
| Electrodes: See LEEP | |
| Endocervical Samplers | |
| Cytobrush | 65 |
| Endocervical Specula | |
| LEEP | 27 |
| Stainless Steel | 46-47 |
| Foerster Sponge Forceps | 47 |
| Endometrial Samplers | |
| Disposable..... | 66 |
| Endosee® Advance | 6-7 |

| | |
|--|-------|
| Euro-Med Instruments (see also LEEP) | |
| Biopsy Punches | 38-41 |
| Needle Holders | 58 |
| Sounds and Dilators | 50-51 |
| Stainless Steel Curettes..... | 52-53 |
| Stainless Steel Endocervical Specula | 46-47 |
| Stainless Steel Forceps..... | 56-58 |
| Stainless Steel Hooks..... | 54 |
| Stainless Steel Lateral Wall Retractors..... | 48 |
| Stainless Steel Scissors | 59 |
| Stainless Steel Specula..... | 42-45 |
| Stainless Steel Tenacula | 55 |
| Fertility and Family Planning | |
| Fertility Assessment | 77 |
| CBR® Newborn Stem Cell Preservation | 80 |
| Aspirette Endocervical Aspirator | 77 |
| H/S Catheter..... | 76 |
| H/S Elliptosphere Catheter | 76 |
| Fetal Monitors | |
| F9..... | 81 |
| Fetal 2EMR | 81 |
| Forceps | |
| LEEP | 26 |
| Stainless Steel | 56-58 |
| General Examination: Diagnostic Disposables | |
| Breast Aspiration Dual-Port Biopsy Needle..... | 68 |
| Catheters | 69 |
| Cervical and Endocervical Samplers..... | 65 |
| Dilators..... | 64 |
| Milex Cervical Spatulas..... | 69 |
| Word Catheter..... | 69 |
| Gynecological Specialty Instruments | |
| Dilators..... | 50-51 |
| Endometrial Curettes | 53 |
| Euro-Med Stainless Steel Endocervical Curettes..... | 52 |
| Euro-Med Stainless Steel Endocervical Specula..... | 46-47 |
| Euro-Med Stainless Steel Hooks..... | 54 |
| Euro-Med Stainless Steel Specula | 42-45 |
| Stainless Steel Forceps..... | 56-58 |
| Lateral Wall Retractors | 48 |
| Uterine Sounds..... | 50 |

| | |
|--|-------|
| H/S Elliptosphere Catheter | 76 |
| Hooks | |
| LEEP | 26-27 |
| Stainless Steel..... | 54 |
| Hysteroscopy | |
| Endosee® Advance..... | 6-7 |
| Lateral Wall Retractors | |
| LEEP | 24 |
| Stainless Steel..... | 48 |
| LEEP | |
| Accessories..... | 20-21 |
| Electrodes | 16-17 |
| Endocervical Specula | 27 |
| Fischer Cone Biopsy Excisor | 18 |
| Forceps..... | 26 |
| Hooks | 26-27 |
| Lateral Wall Retractors..... | 24 |
| LEEP PRECISION™ Integrated System | 14-15 |
| Procedure Essentials..... | 19 |
| Retractors | 24 |
| Rulers..... | 27 |
| Tenacula..... | 26 |
| Tubing | 20-21 |
| Vaginal Specula..... | 22-23 |
| Wallach® Quantum 2000 Accessories..... | 28-29 |
| LEISEGANG | |
| Colposcopy | 30-33 |
| LM-900 Cryosurgery System..... | 70-71 |
| Lighted Specula and Other Lighted Devices | |
| ER-Spec™ | 36 |
| Her-Spec® | 36 |
| OfficeSpec® | 37 |
| Anospec® | 37 |
| Milex | |
| Dilators..... | 64 |
| Pessaries..... | 10-13 |
| Needle Holders | 58 |
| Needles | |
| Potocky Needle | 19 |
| Obstetrical Care | 81-83 |
| Ophthalmic Cryosurgery | 90 |

| | |
|--|------------|
| Pessaries: See Urinary Incontinence | |
| Paragard | 78-79 |
| Pipelle | 60 |
| Pipet Curet | 60 |
| Retractors | |
| LEEP | 24 |
| Stainless Steel | 43, 48 |
| Scissors – Stainless | 59 |
| Sounds | 50, 63, 67 |
| Specula | |
| Dimensions | 44-45 |
| LEEP | 22-23 |
| Stainless Steel | 42-43 |
| Tenacula | |
| LEEP | 26 |
| Stainless Steel | 55 |
| Tubal Patency | |
| ABBI..... | 74 |
| Tubing | |
| Smoke Evacuator Tubing | 20 |
| Tubing Set | 21 |
| Urinary Incontinence | |
| Milex Pessaries | 10-13 |
| Trimo-San | 10 |
| Uterine Infusion Catheters | |
| H/S Catheters..... | 76 |
| H/S Elliptosphere Catheter | 76 |
| Uterine Sounds | 50 |
| Uterine Structure Evaluation | |
| ABBI | 74 |
| Water Vapor Ablation | |
| Mara | 8-9 |

Committed to Our Customers

At CooperSurgical, we strive for the highest level of expertise and knowledge of our products and to develop an ongoing relationship with each of our customers that is built on trust, quality and consistency. We look to provide service that exceeds not only our customers’ expectations, but industry standards as well. CooperSurgical’s Customer Service Department provides full-service order entry with extended hours to facilitate order placement. As an added convenience we also now offer online order entry. Call now to have your Customer Service Representative set up your account today. Once you are registered for online ordering, you will be able to place orders, track shipments, review prior purchase history and schedule orders for future delivery.

- Available for sale where CE mark is required
- Available for sale in Canada

CE and CA labeling is subject to change. Contact Customer Service or your local sales representative for the most up-to-date information.

Warranty Information

CooperSurgical, Inc. ("CooperSurgical") warrants that its products offered in this catalog, **other than** single-use products and consumables, will be free from defects in materials and workmanship for one year from the date of purchase. CooperSurgical warrants that its single-use products and consumables offered in this catalog will, through the earlier of 90 days from the date of purchase of each such product and the expiration date of that product stated on that product's package, be free from defects in materials and workmanship that are discovered prior to or during the single use of that product.

If CooperSurgical determines that a product offered in this catalog fails to conform to the CooperSurgical warranty applicable to that product, CooperSurgical will, as the sole remedy for that failure to conform, repair or replace that product, at CooperSurgical's discretion, free of charge.

To return a product to CooperSurgical, a customer must comply with CooperSurgical's Returned Goods Policy described in this catalog. A customer will not have any remedy if a CooperSurgical product does not conform to the warranty for that product unless that product is returned to CooperSurgical in accordance with that Returned Goods Policy. CooperSurgical will ship returned products that it repaired or replaced under warranty to the customer who returned those products, at CooperSurgical's expense F.O.B. the customer's facility. Under all other circumstances, CooperSurgical will ship returned products to the customer who returned those products at the customer's expense F.O.B. CooperSurgical's facility.

CooperSurgical's warranties, including those in this catalog, do not cover damage caused by misuse, improper care, improper use of chemicals or cleaning methods, loss, theft, servicing or sharpening by non-authorized personnel or negligent or intentional conduct on the part of the owner or user of the product. CooperSurgical's warranties do not cover normal wear and tear or general maintenance, including cleaning and sharpening services required as a result of normal use. Any modifications or changes to a product will void that product's warranty.

THE ABOVE WARRANTIES ARE IN LIEU OF, AND COOPERSURGICAL HEREBY DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, WRITTEN OR ORAL, WITH RESPECT TO COOPERSURGICAL'S PRODUCTS, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO TERMS, CONDITIONS, UNDERSTANDINGS OR AGREEMENTS THAT PURPORT TO MODIFY THE ABOVE WARRANTIES OR THAT MAKE ANY ADDITIONAL WARRANTIES FOR ANY COOPERSURGICAL PRODUCT SHALL HAVE ANY LEGAL EFFECT UNLESS MADE IN WRITING AND SIGNED BY AN AUTHORIZED COOPERSURGICAL CORPORATE OFFICER.

COOPERSURGICAL SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR LOST PROFITS, DAMAGES FROM LOSS OF USE OR LOST DATA, OR INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ITS WARRANTIES OR OTHERWISE FOR ANY CLAIM RELATED TO COOPERSURGICAL'S PRODUCTS, EVEN IF COOPERSURGICAL HAS BEEN ADVISED, KNEW OR SHOULD HAVE KNOWN OF THE POSSIBILITY OF SUCH DAMAGES. COOPERSURGICAL'S LIABILITY WITH RESPECT TO A PRODUCT COVERED BY A WARRANTY OR OTHERWISE SHALL BE LIMITED IN ALL CIRCUMSTANCES TO THE PURCHASE PRICE OF THAT PRODUCT.

5 Ways to Order

- Local Sales Representative**
www.coopersurgical.com/pages/salesrepsearch.aspx
- Email**
U.S.: Orders@coopersurgical.com
Int'l: InternationalCustomerservice@coopersurgical.com
- E-Ordering (U.S. Only)**
Call 800.243.2974 to set up your account today.
- Phone Orders**
800.243.2974 (8:00 AM – 8:30 PM EST)
Int'l: 203.601.9818
- Fax Orders**
800.262.0105 (24 hours a day)
Int'l: 203.601.4747

Returned Goods Policy

- Goods will be accepted for return for the following reasons:**
 - If shipment was made without Customer's authorization or order
 - If incorrect goods were shipped
 - If defective goods are covered by the standard warranty
- To return goods, you must contact a Customer Service Representative by telephone (1.800.243.2974 or 1.203.601.5200) for a Returned Merchandise Authorization (RMA) number. Items will not be accepted without an RMA number.
Please have the following information:
 - Reason you wish to return the goods
 - Quantity, description, part number, serial number of the goods
 - Date of Receipt of Order
 - Customer's Purchase Order and the CooperSurgical Invoice Number**All used instruments must be cleaned and sterilized prior to being shipped for return.**
- Shipment must be sent prepaid by Customer. Freight collect shipments will not be accepted, and goods will be returned to sender.
- If Customer intends to return equipment ordered in error, the following restocking charges and terms will apply:
 - There is no restocking fee if goods are returned within 30 days of shipment
 - There is a restocking fee of 25% if goods are returned before 31-60 days of shipment
 - Goods must be returned unused, in the original carton, and in marketable condition
 - Refurbishing and replacement charges will be added to the restocking charges for damaged or missing items
 - No return after 60 days
 - No refund on sterile, single-use disposable products

Send All Returns to: **CooperSurgical**
95 Corporate Drive
Trumbull, CT 06611

Trademark Notices

All trademarks are the property of their respective owners.

Endosee® Advance

Mara®

Milex®

LEEP PRECISION™

Fischer Cone Biopsy Excisor®

Pipelle®

Leisegang™

ER-Spec®

Her-Spec®

Euro-Med™

H/S Elliptosphere®

ABBI® Air Bubble Based Infuser

Paragard® Intrauterine
Copper Contraceptive

Cord Blood RegistrySM

Wallach®



CooperSurgical®

800.243.2974

www.coopersurgical.com

95 Corporate Drive, Trumbull, CT 06611