

Her-Spec[®]

Single-Use Click-Free Vaginal Speculum with Integrated LED Light Source



Disposable, Bright Built-In Light, Batteries Pre-Installed

Her-Spec utilizes a CLICK-FREE design that eliminates the clicking sound associated with the ratchet mechanism in other disposable vaginal specula.

Ideal for clinics, emergency rooms, urgent care centers, and physician offices.

Her-Spec is designed for use on women undergoing a procedure requiring vaginal access and exposure.

Product Features

- Disposable, single-use device reduces risk of cross-contamination
- Each Her-Spec comes ready for use with an integrated bright cordless LED light source
- Single-use, battery-operated LED light source can produce peak light for 20+ minutes
- Smooth locking feature allows clinicians to set blades to desired open position by pushing center tab to automatically lock blades in place
- Ergonomically designed handle with diamond grip for added comfort and maneuverability during exams
- 105-degree angled handle creates extra space between clinician's hand and patients' body and/or exam table
- Touch free battery removal feature allows users the option to easily eject batteries without contact after use



CooperSurgical[®]

Her-Spec[®]

Disposable, Bright Built-In Light, Batteries Pre-Installed

CooperSurgical's lighted vaginal speculum comes with pre-installed batteries, and offer a wide variety of styles and sizes to best meet the needs of your practice.



Ordering Information

NUMBER	DESCRIPTION
C020500	Her-Spec [®] (Small)
C020510	Her-Spec [®] (Medium)
C020520	Her-Spec [®] (Large)
C020600	Her-Spec [®] (Small - Sterile)
C020610	Her-Spec [®] (Medium - Sterile)
C020620	Her-Spec [®] (Large - Sterile)



800.243.2974 | 203.601.5200
www.coopersurgical.com



CooperSurgical[®]

© 2021 CooperSurgical, Inc.

C-US-OBP-000007